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Pharmacopeial Convention, Inc., 12601 Twinbrook Parkway, Rockville, MD

Fig. 5

## Instructions for Blood Collection (CPDA-1 or CPD/ADSOL™) BLOOD-PACK™ Unit

### Rx only

Before proceeding further, identify BLOOD-PACK™ unit using appropriate donor identification system.

1. Suspend the BLOOD-PACK unit as far below donor's arm as possible.
2. Apply blood pressure cuff or tourniquet.
3. Disinfect site of venipuncture.
4. Tie tourniquet or inflate pressure cuff to 60 mm Hg.
5. Appropriately secure donor needle and/or tubing.
6. **Clamp donor tube**, remove needle cover and perform venipuncture.
7. Release clamp from donor tube to permit flow of blood.
8. **Mix blood and anticoagulant at several intervals during collection and immediately after collection.**
9. Collect the quantity of blood stated on the pack label. Anticoagulant is sufficient for  $\pm 10\%$ .
10. Hermetically seal the tube between the sampling site and primary collection bag approx. 3 cm from sample site. (A heat sealer or metal clip or other appropriate method may be used to seal at this point). Strip the donor tubing  $\pm 10$  cm from the clamp toward the collection pack, bend the tubing and apply a metal clip (4R4418) onto the stripped portion of the tubing. To collect samples assemble a Luer needle adapter and tube holder.
11. Insert the assembly clockwise firmly and positively into the sample site ensuring the Luer is fully entered into the sample site. Proceed to blood sampling with evacuated sample tubes.
12. Isolate the BLOOD-PACK by cutting the tubing (or breaking the seal) between the hemostat and metal clip.
13. Remove tourniquet or release pressure in cuff and withdraw needle. Dispose of needle assembly in an appropriate biohazardous waste container following established procedures.
14. If necessary, strip blood from donor tube into the BLOOD-PACK unit, mix and allow tube to refill. Repeat once. Seal at X marks, beginning at sealed end of donor tube toward pack, to provide numbered aliquots of anticoagulated blood for crossmatching.
15. Remove and discard the Y-Sampling Site and the DONORCARE™ Needle Guard into an appropriate biohazardous waste container following established procedures.
16. **Component Preparation (Manual separation):**
  - Fresh frozen plasma should be separated from the red blood cells and placed in the freezer at  $-18^{\circ}\text{C}$  or colder within 8 hours after blood collection.

- If applicable, ADSOL™ red cell preservation solution should be added to the red blood cells immediately after the removal of plasma. Preparation of AS-1 red blood cells may vary depending on processing option selected:
    - a) Within 8 hours of blood collection if whole blood is held at ambient temperature.
    - b) Within 3 days of blood collection if whole blood is refrigerated.
  - For component preparation using manual plasma extractors, use standard component processing and storage techniques.
  - If a platelet concentrate is to be prepared, it should be separated from the red blood cells within 8 hours after blood collection.
17. When processing a multiple BLOOD-PACK unit, centrifuge primary and secondary containers to prepare CPD red blood cells, or CPDA-1 red blood cells as applicable.
  18. Place primary container in plasma extractor and express plasma into empty TRANSFER PACK™ container by releasing pressure plate and opening closure in tubing of primary container. (Break the cannula on top of the primary pack by holding the rigid part and **bending 3 times laterally 90°**. Make sure the plasma flows to the TRANSFER PACK™).
  19. When desired amount of plasma has been removed, clamp tubing between Y and plasma container.
  20. If applicable, suspend ADSOL red cell preservation solution container, open closure in tubing (Break the stem of the cannula using the other hand by **bending back and forth at 90° at least 3 times**. Ensure the cannula is well broken by locating the broken part in the upper part of the housing) and drain contents into primary container of CPD red blood cells. Clamp tubing.
  21. Seal transfer tubing in three places after second segment number near primary container (leaving two segment numbers connected to the primary container) and cut middle seal being careful to avoid fluid splatter. **For Double BLOOD-PACK unit codes, discard ADSOL solution container. For other ADSOL codes**, the empty solution container may be used as a TRANSFER PACK container for further component preparation.
  22. If applicable, mix ADSOL red cell preservation solution and red cells thoroughly.
  23. **Component Preparation (OPTIPAC™/OPTIPRESS™):**
    - Refer to OPTIPRESS™ manual for detailed instructions.
    - For component preparation using an OPTIPAC™, use standard component processing and storage techniques.
  24. Centrifuge the pack to separate blood component.
  25. Position the OPTIPAC primary pack in the OPTIPRESS.
  26. Put the additive solution container on the bench and the tubing through the lower clamp.
  27. Put the empty TRANSFER PACK containers on top of the OPTIPRESS and tubing through the upper clamp.
  28. Start the OPTIPRESS by pushing the “Press” button.
  29. Break the cannula on top of the primary pack by holding the rigid part and **bending 3 times laterally 90°**. Make sure the plasma flows to the TRANSFER PACK.
  30. As soon as the top clamp is closed, hold the base of the cannula on the additive solution container between the thumb and first finger. Break the stem of the cannula using the other hand by **bending back and forth at 90° at least 3 times**. Ensure the cannula is well broken by locating the broken part in the upper part of the housing.
  31. When the procedure is completed, both the upper and lower clamps occlude the tubing so as to prevent red cell contamination. At this point place a clip on the line at the top of the OPTIPAC to prevent any red cells travelling up the plasma line when the OPTIPAC has been removed from the

OPTIPRESS.

32. Seal the TRANSFER PACK containing plasma and the unit of red cells, leaving behind the OPTIPAC and the platelet pack.
33. Push the “PRESS” button on the OPTIPRESS to release the pressure plate and remove the packs.
34. After storage of the OPTIPAC, the primary pack, and the empty platelet container may be used to prepare other components according to blood center procedures.
35. **Storage:**
  - Store suspended AS-1 red blood cells or suspended CPDA-1 whole blood/red blood cells (as applicable) between 1 and 6°C.
  - Infuse AS-1 red blood cells within 42 days of collection, and CPDA-1 whole blood/red blood cells within 35 days of collection.

#### NOTES:

- Protect container and tubing from sharp objects.
- When frozen, plastic is more fragile.

#### CAUTIONS TO BE TAKEN

- Do not use if there is visible sign of deterioration.
- Do not use if fluid path closures are loose or not intact.
- Do not use unless solution is clear and no leaks are detected.

Sterile, non-pyrogenic fluid path. Steam sterilized.

Single use only. Do not vent.

Dispose of container appropriately.

Store at Controlled Room Temperature. USP Definition of “Controlled Room Temperature” United States Pharmacopeia, General Notices.

United States Pharmacopeial Convention, Inc.

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Imported and distributed in Indonesia by:

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Imported and distributed in Thailand by:

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Imported and distributed in Perú by:

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Imported and distributed in Venezuela by:

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Imported and distributed in Colombia by:

**Fenwal Colombia Ltda.**

Centroempresa, Cali, Colombia

Reg. No.: INVIMA-2006M-0000140-R1

07-19-04-576 REV: A 03/2010

FENWAL, BOLSANG, ADSOL, TRANSFER PACK, OPTIPRESS y OPTIPAC son Marcas Registradas de Fenwal, Inc.

FENWAL, BLOOD-PACK, ADSOL, TRANSFER PACK, OPTIPRESS and OPTIPAC are trademarks of Fenwal, Inc.

DONORCARE is a trademark of ITL Corporation.

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**Plastic**

**Rx Only**

Quadruple BLOOD-PACK™ OPTIPAC™ unit consisting of a primary pack containing 63 mL of CPD anticoagulant solution for collection of 450 mL of blood, one TRANSFER PACK™ containing 100 mL of ADSOL™ solution and two TRANSFER PACK™ containers without solution. 16 gauge needle.

See instructions for use. Sterile, non-pyrogenic fluid path. Steam sterilized. Single use only. Do not vent. Do not use if there is any visible evidence of deterioration. Dispose of container appropriately.

Store at Controlled Room Temperature (refer to direction insert).

Unused packs in open foil pouch may be kept 60 days by folding and securing open end of foil pouch, to prevent possible loss of moisture.

Direct handling of product surfaces prior to extended storage in the foil pouch may result in mold growth.

Units removed from the foil pouch must be used within 4 days (96 hours). Units out of the foil pouch for longer than 4 days must be discarded.

Each 63 mL of CPD solution contains: 1.66 g, Sodium citrate (dihydrate), USP, 1.61 g, Dextrose (monohydrate), USP, 188 mg, Citric acid (anhydrous), USP, 140 mg Monobasic sodium phosphate (monohydrate), USP. Water for Injection USP to 63 mL.

Each 100 mL of ADSOL solution contains: Sodium chloride, USP, 900 mg, Dextrose monohydrate, USP, 2.20 g, Adenine, USP, 27.0 mg, Mannitol, USP, 750 mg. Water for Injection USP to 100 mL.

Manufactured by:

**Fenwal International, Inc.**

Road 357, Km. 0.8

Maricao, PR 00606

Made in USA

FENWAL, ADSOL, OPTIPAC, TRANSFER PACK and BLOOD-PACK are trademarks of Fenwal, Inc.

07-28-04-571 REV: A

Imported and distributed in India by:

**Fenwal India Pvt Ltd**

Upper Ground Floor, Tower B

DLF Building No. 10

DLF Cyber City, DLF Phase-II, Gurgaon 122 002

Haryana, India

Import License No.: FF-504-14890

MRP Rs.: 850 (Inclusive of all taxes)

Lot No.:

Exp.:

R4R7794

4 Units

# QUADRUPLE BLOOD-PACK™ OPTIPAC™ UNIT, CPD/ADSOL™

PL 146  
PL 1240 Plastic

Rx Only

Quadruple BLOOD-PACK™ OPTIPAC™ unit consisting of a primary pack containing 63 mL of CPD anticoagulant solution for collection of 450 mL of blood, one TRANSFER PACK™ containing 100 mL of ADSOL™ solution and two TRANSFER PACK™ containers without solution. 16 gauge needle.

See instructions for use. Sterile, non-pyrogenic fluid path. Steam sterilized. Single use only. Do not vent. Do not use if there is any visible evidence of deterioration. Dispose of container appropriately.

Store at Controlled Room Temperature (refer to direction insert).

Unused packs in open foil pouch may be kept 60 days by folding and securing open end of foil pouch, to prevent possible loss of moisture.

Direct handling of product surfaces prior to extended storage in the foil pouch may result in mold growth.

Units removed from the foil pouch must be used within 4 days (96 hours). Units out of the foil pouch for longer than 4 days must be discarded.

Each 63 mL of CPD solution contains: 1.66 g, Sodium citrate (dihydrate), USP, 1.61 g, Dextrose (monohydrate), USP, 188 mg, Citric acid (anhydrous), USP, 140 mg Monobasic sodium phosphate (monohydrate), USP. Water for Injection USP to 63 mL.

Each 100 mL of ADSOL solution contains: Sodium chloride, USP, 900 mg, Dextrose monohydrate, USP, 2.20 g, Adenine, USP, 27.0 mg, Mannitol, USP, 750 mg. Water for Injection USP to 100 mL.

 Manufactured by:  
**Fenwal International, Inc.**  
Road 357, Km. 0.8  
Maricao, PR 00606

Made in USA

FENWAL, ADSOL, OPTIPAC, TRANSFER PACK and BLOOD-PACK are trademarks of Fenwal, Inc.

07-28-04-571 REV: A

Imported and distributed in India by:

**Fenwal India Pvt Ltd**

Upper Ground Floor, Tower B

DLF Building No. 10

DLF Cyber City, DLF Phase-II, Gurgaon 122 002

Haryana, India

Import License No.: FF-504-14890

MRP Rs.: 850 (Inclusive of all taxes)

Lot No.:

Exp.:

## ADSOL RED CELL PRESERVATION SOLUTION SYSTEM IN PLASTIC CONTAINER (PL 146 PLASTIC)

anticoagulant citrate phosphate dextrose (cpd) solution and adsol preservation solution kit

### Product Information

<b>Product Type</b>	HUMAN PRESCRIPTION DRUG	<b>Item Code (Source)</b>	NDC:0942-6494
---------------------	-------------------------	---------------------------	---------------

### Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:0942-6494-04	1 in 1 KIT		

### Quantity of Parts

Part #	Package Quantity	Total Product Quantity
Part 1	1 BAG	63 mL

Part 2 | 1 BAG

100 mL

**Part 1 of 2****CPD**

citrate phosphate dextrose solution

**Product Information**

<b>Route of Administration</b>	INTRAVENOUS
--------------------------------	-------------

**Active Ingredient/Active Moiety**

Ingredient Name	Basis of Strength	Strength
<b>Trisodium Citrate Dihydrate</b> (UNII: B22547B95K) (Anhydrous Citric Acid - UNII:XF417D3PSL)	Anhydrous Citric Acid	1.66 g in 63 mL
<b>Dextrose Monohydrate</b> (UNII: LX22YL083G) (ANHYDROUS DEXTROSE - UNII:5SL0G7R0OK)	Dextrose Monohydrate	1.61 g in 63 mL
<b>Anhydrous Citric Acid</b> (UNII: XF417D3PSL) (Anhydrous Citric Acid - UNII:XF417D3PSL)	Anhydrous Citric Acid	188 mg in 63 mL
<b>Sodium Phosphate, Monobasic, Monohydrate</b> (UNII: 593YOG76RN) (PHOSPHATE ION - UNII:NK08V8K8HR)	Sodium Phosphate, Monobasic, Monohydrate	140 mg in 63 mL

**Inactive Ingredients**

Ingredient Name	Strength
<b>Sodium Hydroxide</b> (UNII: 55X04QC32I)	
<b>Water</b> (UNII: 059QF0KO0R)	

**Packaging**

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1		63 mL in 1 BAG; Type 0: Not a Combination Product		

**Marketing Information**

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
NDA	BN8 11104	05/25/2010	

**Part 2 of 2****ADSOL RED CELL PRESERVATION SOLUTION SYSTEM**

adsol red cell preservation solution solution

## Product Information

Route of Administration INTRAVENOUS

## Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
<b>Dextrose Monohydrate</b> (UNII: LX22YL083G) (ANHYDROUS DEXTROSE - UNII:5SL0G7R0OK)	Dextrose Monohydrate	2.2 g in 100 mL
<b>Sodium Chloride</b> (UNII: 451W47IQ8 X) (SODIUM CATION - UNII:LYR4M0NH37)	Sodium Chloride	900 mg in 100 mL
<b>Mannitol</b> (UNII: 3OWL53L36A) (Mannitol - UNII:3OWL53L36A)	Mannitol	750 mg in 100 mL
<b>Adenine</b> (UNII: JAC85A2161) (Adenine - UNII:JAC85A2161)	Adenine	27 mg in 100 mL

## Inactive Ingredients

Ingredient Name	Strength
<b>Water</b> (UNII: 059QF0K00R)	

## Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1		100 mL in 1 BAG; Type 0: Not a Combination Product		

## Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
NDA	BN8 11104	05/25/2010	

## Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
NDA	BN8 11104	05/25/2010	

**Labeler** - Fenwal, Inc. (794519020)

## Establishment

Name	Address	ID/FEI	Business Operations
Fenwal International, Inc.		091164590	MANUFACTURE(0942-6494)

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© 2010 Fenwal, Inc. Todos los derechos reservados. All rights reserved.