

**TERUFLEX BLOOD BAG SYSTEM ANTICOAGULANT CITRATE PHOSPHATE
DEXTROSE (CPD) AND OPTISOL RED CELL PRESERVATIVE - anticoagulant citrate
phosphate dextrose (cpd) and as-5 red cell preservative
Terumo Corporation**

**TERUFLEX® BLOOD BAG SYSTEM CPD/OPTISOL® SOLUTION BLOOD BAG SYSTEM
TERUFLEX® BLOOD BAG SYSTEM
CPD/OPTISOL® SOLUTION BLOOD BAG SYSTEM**

Read these instructions carefully before use. Rx ONLY.

INTENDED USE

Collection of 250 mL Whole Blood:

- From elderly, pediatric, or smaller autologous donors who may not tolerate larger blood collections.
- For blood collections as directed by a patient's physician.
- Platelets prepared from 250 mL whole blood collections can only be used for the autologous donation patient.
- Platelet preparation from a 250 mL whole blood collection are unlikely to meet the 5.5×10^{10} platelet concentration required by the regulations.

INSTRUCTIONS FOR BLOOD COLLECTION: Use aseptic technique

1. Confirm that all numbered tubing of each blood bag unit contains its own identical segment numbers.
2. Make a loose knot in the donor tubing approximately 10 cm from the needle.
3. Clamp donor tubing with hemostat.
4. Suspend primary bag a minimum of 60 cm below the donor's arm
5. Apply blood pressure cuff. Disinfect site of phlebotomy. Inflate blood pressure cuff to 60 mmHg.
6. Remove needle protector and perform phlebotomy. Remove hemostat to permit blood flow into primary bag.

CAUTION Do not touch needle after removing the needle protector.

7. Tape donor tubing securely to donor's arm.
8. MIX BLOOD WITH ANTICOAGULANT AT SEVERAL INTERVALS DURING COLLECTION.
9. Collect 250 mL of blood.

10. Tighten knot firmly after collection. Clamp between knot and needle. Sever donor tubing between knot and clamp. Collect pilot samples.

11. Reapply clamp to donor tubing; release blood pressure cuff and remove needle from donor's arm.

CAUTION Discard tubing/phlebotomy needle unit according to institutional procedures.

12. Immediately after collection, invert bag several times mixing blood with anticoagulant thoroughly.
13. Strip blood from donor tubing into bag, mix well, and allow tubing to refill. Seal at X marks on donor tubing to provide numbered aliquots of anticoagulated blood for testing.

CAUTION Begin sealing at needle end and work towards bag.

14. The time of addition of OPTISOL solution may vary depending on the processing option selected. Add solution under one of the following conditions.

- a) After removal of plasma from freshly collected blood.
- b) Within 8 hours of blood collection if platelets are prepared.
- c) Within 72 hours of collection if blood is refrigerated immediately following collection.

15. Centrifuge unit to separate red cells from plasma.

16. Snap CLIKTIP (incline closure device) of primary collection bag and transfer plasma into satellite bag. Clamp transfer tubing of satellite bag.

17. Snap CLIKTIP of OPTISOL solution bag and drain contents into primary bag containing red blood cells. Seal tubing of primary bag in two places, cut between seals, and separate from satellite bag.

NOTE: For TERUFLEX double bags, seal OPTISOL solution bag tubing in two places, but between

seals, and separate. Discard OPTISOL solution container.


18. Mix OPTISOL solution and red cells thoroughly. Store between 1-6°C.

19. Infuse OPTISOL red blood cells within 42 days of collection.

For further processing, use standard component processing techniques.

To open the blister package, peel the cover film back 4/5 of its length. (See Fig. 3)

After opening, unused blood bags may be stored at room temperature for 96 hours, or they may be stored for 30 days in the blister package after returning the cover film to the original position and sealing with tape to prevent possible loss of moisture.



CAUTIONS

- THE PACKET OF AGELESS CONTAINED IN THIS PACKAGE GENERATES HEAT UPON REMOVAL AND SHOULD BE HANDLED WITH CARE.
 - DISPOSE WITH AGELESS PACKET IN TRAY.
 - DO NOT DISPOSE WITH WASTES CONTAINING VOLATILE OR FLAMMABLE MATERIALS.
 - DISCARD AGELESS PACKET WITHOUT OPENING.
- DO NOT USE UNLESS ANTICOAGULANT IS CLEAR.
RECOMMENDED STORAGE: Room Temperature (25°C/77°F).
Avoid excessive heat and direct sunlight. Protect from freezing.

TERUMO CORPORATION

44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO, JAPAN

TERUMO and TERUFLX are registered in U.S. Patent and Trademark Office by TERUMO CORPORATION.

OPTISOL and CLIKTIP are trademarks of TERUMO CORPORATION.

Issued 2/95

N-BB-OP-A3

Tray/Case Label

TERUFLEX® BLOOD BAG SYSTEM with Blood Sampling Arm®

CPD WITH OPTISOL® RED CELL PRESERVATIVE SOLUTION FOR COLLECTION OF 250 mL OF BLOOD

Each unit consists of a primary bag containing 35mL of Anticoagulant CPD solution, with a satellite bag containing 56mL of OPTISOL Red Cell Preservative Solution.

Each 35 mL Anticoagulant CPD solution USP contains 893mg Dextrose (monohydrate) USP, 921mg Sodium Citrate (dihydrate) USP, 105mg Citric Acid (anhydrous) USP, 77.7mg Monobasic Sodium Phosphate (monohydrate) USP.

Each 56 mL OPTISOL Red Cell Preservative Solution contains 491mg Sodium Chloride USP, 504mg Dextrose (monohydrate) USP, 294mg Mannitol USP, 16,8mg Adenine USP.

STERILE, NON-PYROGENIC FLUID PATH.
DO NOT USE UNLESS ANTICOAGULANT IS CLEAR

CODE

LOT No.

EXPIRY

UNITS

DONOR NEEDLE 17G x 1 1/2" (1.40 x 38mm)

Rx ONLY

RECOMMENDED STORAGE: Room Temperature (15-30°C/59-86°F).

Avoid excessive heat. Protect from freezing.

After opening, unused bags may be stored for 30 days by returning cover film to original position and sealing with tape to prevent possible loss of moisture. See Instructions For Blood Collection.

Manufactured by : **TERUMO CORPORATION** Tokyo, Japan

® : Registered Trademark of TERUMO CORPORATION



TERUFLEX[®] BLOOD BAG SYSTEM with Blood Sampling Arm[®]

**CPD WITH OPTISOL[®] RED CELL PRESERVATIVE SOLUTION
FOR COLLECTION OF 250mL OF BLOOD**

Each unit consists of a primary bag containing 35mL of Anticoagulant CPD solution, with a satellite bag containing 58mL of OPTISOL Red Cell Preservative Solution.

Each 35mL Anticoagulant CPD solution USP contains 893mg Dextrose (monohydrate) USP, 921mg Sodium Citrate (dihydrate) USP, 105mg Citric Acid (anhydrous) USP, 77.7mg Monobasic Sodium Phosphate (monohydrate) USP.

Each 58mL OPTISOL Red Cell Preservative Solution contains 491mg Sodium Chloride USP, 504mg Dextrose (monohydrate) USP, 294mg Mannitol USP, 16.8mg Adenine USP.

STERILE, NON-PYROGENIC FLUID PATH.

DO NOT USE UNLESS ANTICOAGULANT IS CLEAR.

CODE

LOT No.

EXPIRY

UNITS

DONOR NEEDLE 17G x 1½" (1.40 x 38mm)

Rx ONLY

RECOMMENDED STORAGE : Room Temperature (15-30 °C / 59-86 °F).

Avoid excessive heat. Protect from freezing.

After opening, unused bags may be stored for 30 days by returning cover film to original position and sealing with tape to prevent possible loss of moisture.

See Instructions For Blood Collection.

Manufactured by : TERUMO CORPORATION Tokyo, Japan

® : Registered Trademark of TERUMO CORPORATION

Rev. 07/04

B-4-B7-A4 ①

**TERUFLEX BLOOD BAG SYSTEM ANTICOAGULANT CITRATE PHOSPHATE
DEXTROSE (CPD) AND OPTISOL RED CELL PRESERVATIVE**

anticoagulant citrate phosphate dextrose (cpd) and as-5 red cell preservative kit

Product Information

Product Type	HUMAN PRESCRIPTION DRUG	Item Code (Source)	NDC:53877-005
---------------------	-------------------------	---------------------------	---------------

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:53877-005-01	30 in 1 CASE		
1		1 in 1 BAG; Type 9: Other Type of Part 3 Combination Product (e.g., Drug/Device/Biological Product)		

Quantity of Parts

Part #	Package Quantity	Total Product Quantity
Part 1	1 BAG	35 mL
Part 2	1 BAG	56 mL

Part 1 of 2**ANTICOAGULANT CITRATE PHOSPHATE DEXTROSE (CPD)**

anticoagulant citrate phosphate dextrose (cpd) solution

Product Information

Route of Administration	INTRAVENOUS
--------------------------------	-------------

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
Trisodium Citrate Dihydrate (UNII: B22547B95K) (ANHYDROUS CITRIC ACID - UNII:XF417D3PSL)	ANHYDROUS CITRIC ACID	26.3 g in 1000 mL
SODIUM PHOSPHATE, MONOBASIC, UNSPECIFIED FORM (UNII: 3980JIH2SW) (PHOSPHATE ION - UNII:NK08V8K8HR, SODIUM CATION - UNII:LYR4M0NH37)	SODIUM PHOSPHATE, MONOBASIC, UNSPECIFIED FORM	2.22 g in 1000 mL
Dextrose Monohydrate (UNII: LX22YL083G) (ANHYDROUS DEXTROSE - UNII:5SL0G7R0OK)	Dextrose Monohydrate	25.5 g in 1000 mL
Anhydrous Citric Acid (UNII: XF417D3PSL) (Anhydrous Citric Acid - UNII:XF417D3PSL)	Anhydrous Citric Acid	2.99 g in 1000 mL

Inactive Ingredients

Ingredient Name	Strength
Water (UNII: 059QF0K00R)	

Packaging

#	Item	Package Description	Marketing Start	Marketing End
---	------	---------------------	-----------------	---------------

#	Code	Package Description	Date	Date
1		35 mL in 1 BAG; Type 9: Other Type of Part 3 Combination Product (e.g., Drug/Device/Biological Product)		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
NDA	BN880217	12/15/2009	

Part 2 of 2

OPTISOL RED CELL PRESERVATIVE

as-5 red cell preservative solution

Product Information

Route of Administration INTRAVENOUS

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
Sodium Chloride (UNII: 451W47IQ8X) (SODIUM CATION - UNII:LYR4M0NH37)	Sodium Chloride	877 mg in 100 mL
Dextrose Monohydrate (UNII: LX22YL083G) (ANHYDROUS DEXTROSE - UNII:5SL0G7R0OK)	Dextrose Monohydrate	900 mg in 100 mL
Mannitol (UNII: 3OWL53L36A) (Mannitol - UNII:3OWL53L36A)	Mannitol	525 mg in 100 mL
Adenine (UNII: JAC85A2161) (Adenine - UNII:JAC85A2161)	Adenine	30 mg in 100 mL

Inactive Ingredients

Ingredient Name	Strength
Water (UNII: 059QF0K00R)	

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1		56 mL in 1 BAG; Type 9: Other Type of Part 3 Combination Product (e.g., Drug/Device/Biological Product)		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
NDA	BN880217	12/15/2009	

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
NDA	BN880217	12/15/2009	

Labeler - Terumo Corporation (690543319)

Establishment

Name	Address	ID/FEI	Business Operations
Terumo Corp. - Fujinomiya Factory		695214015	manufacture(53877-005) , STERILIZE(53877-005) , ANALYSIS(53877-005) , LABEL(53877-005)

Revised: 12/2018

Terumo Corporation