



January 8, 2024

InkSpace Imaging, Inc.  
% Taras Bouzakine  
Director, Regulatory Affairs  
Veranex, Inc.  
224 Airport Parkway, Suite 250  
San Jose, California 95110

Re: K233444

Trade/Device Name: InkSpace Imaging Small Body Array  
Regulation Number: 21 CFR 892.1000  
Regulation Name: Magnetic Resonance Diagnostic Device  
Regulatory Class: Class II  
Product Code: MOS  
Dated: December 14, 2023  
Received: December 15, 2023

Dear Taras Bouzakine:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

A handwritten signature in black ink, appearing to read 'D. Krainak', is written over a large, light blue, semi-transparent watermark of the letters 'FDA'.

Daniel M. Krainak, Ph.D.  
Assistant Director  
DHT8C: Division of Radiological Imaging  
and Radiation Therapy Devices  
OHT8: Office of Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

Submission Number (if known)

K233444

Device Name

InkSpace Imaging Small Body Array

Indications for Use (Describe)

The InkSpace Imaging Small Body Array is a receive-only coil, used for obtaining diagnostic images of general human anatomy in adult and pediatric patients such as cardiac, spine, shoulder, elbow, knee, foot, and prostate in Siemens 3.0T magnetic resonance imaging systems. These images, when interpreted by a trained physician, yield information that may assist in diagnosis.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) SUMMARY

---

### 510(k) Notification K233444

#### GENERAL INFORMATION [807.92(A)(1)]

**Applicant:**

InkSpace Imaging, Inc.  
5635 West Las Positas Blvd, Suite 403/404  
Pleasanton, CA 94588  
USA  
Phone: 925-425-7410

**Contact Person:**

Taras Bouzakine  
Director, Regulatory Affairs, Veranex, Inc.  
224 Airport Parkway, Suite 250  
San Jose, CA 95110  
USA  
Phone: 603-205-4874  
Email: [taras.bouzakine@veranex.com](mailto:taras.bouzakine@veranex.com)

**Date Prepared:** November 3, 2023

#### DEVICE INFORMATION [807.92(A)(2)]

**Classification:**

Class II per 21 CFR§892.1000

**Product Code:**

MOS

**Trade Name:**

InkSpace Imaging Small Body Array

**Generic/Common Name:**

Coil, Magnetic Resonance, Specialty

#### PREDICATE DEVICE(S) [807.92(A)(3)]

InkSpace Imaging Body Array (K223487)

#### DEVICE DESCRIPTION [807.92(A)(4)]

The InkSpace Imaging Small Body Array is a phased array, receive-only coil intended to work with Siemens 3.0T MRI scanners to provide optimal signal-to-noise ratio and uniformity for high resolution body MR imaging examinations. The coil is comprised of 24 total channels, with two individual flexible pads, each containing 12 individual elements. The elements are optimized to be flexible to conform to patients' anatomies.

## **510(k) SUMMARY**

---

### **INDICATIONS FOR USE [807.92(A)(5)]**

The InkSpace Imaging Small Body Array is a receive-only coil, used for obtaining diagnostic images of general human anatomy in adult and pediatric patients such as cardiac, spine, shoulder, elbow, knee, foot, and prostate in Siemens 3.0T magnetic resonance imaging systems. These images, when interpreted by a trained physician, yield information that may assist in diagnosis.

### **COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICES [807.92(A)(6)]**

The technological characteristics of the modified InkSpace Imaging Small Body Array are substantially equivalent to the predicate device, InkSpace Imaging Body Array (K223487). Table 1 lists the technological characteristics of the predicate and modified devices and provides rationale to support a determination of substantial equivalence. Any differences between the devices do not raise any new or different questions of safety or effectiveness.

## 510(k) SUMMARY

**Table 1: Summary of Technological Characteristics**

Characteristics	Predicate Device: InkSpace Imaging Body Array (K223487)	Subject Device: InkSpace Imaging Small Body Array	Substantial Equivalence Rationale
<b>Classification</b>	Magnetic Resonance Diagnostic Device per 21 CFR§892.1000	Magnetic Resonance Diagnostic Device per 21 CFR§892.1000	N/A (Same).
<b>Product Code</b>	MOS (Coil, Magnetic Resonance, Specialty)	MOS (Coil, Magnetic Resonance, Specialty)	N/A (Same).
<b>Intended Use</b>	<p>Pediatric Population: The InkSpace Imaging Body Array is a receive-only RF coil designed for use with GE 3.0T MRI systems. It is used for obtaining diagnostic images of the body and cardiac regions of pediatric patients in magnetic resonance imaging systems. The nucleus detected is hydrogen. Anatomic Regions: Cardiac/Body for pediatric patients.</p> <p>Adult Population: The InkSpace Imaging Body Array is a receive-only RF coil designed for use with GE 3.0T MRI systems. It is used for obtaining diagnostic images of general human anatomy such as spine, shoulder, elbow, knee, foot, and prostate for adult patients in magnetic resonance imaging systems. The nucleus detected is hydrogen.</p>	<p>The InkSpace Imaging Small Body Array is a receive-only RF coil designed for use with Siemens 3.0T MRI systems. It is used for obtaining diagnostic images of general human anatomy in adult and pediatric patients such as cardiac, spine, shoulder, elbow, knee, foot, and prostate in magnetic resonance imaging systems. The nucleus detected is hydrogen.</p>	<p>Intended use of the modified device maintains the same general purpose and function as the predicate device, differing only by target MRI system.</p>
<b>Indications for Use</b>	<p>Pediatric Population: The InkSpace Imaging Body Array is a receive-only coil, used for obtaining diagnostic images of pediatric cardiac and body in GE 3.0T magnetic resonance imaging systems. These images, when interpreted by a trained physician, yield information that may assist in diagnosis.</p> <p>Adult Population: The InkSpace Imaging Body Array is a receive-only coil, used for obtaining diagnostic images of general human anatomy in adults, such as</p>	<p>The InkSpace Imaging Small Body Array is a receive-only coil, used for obtaining diagnostic images of general human anatomy in adult and pediatric patients, such as cardiac, spine, shoulder, elbow, knee, foot, and prostate in Siemens 3.0T magnetic resonance imaging systems. These images, when interpreted by a trained physician, yield information that may assist in diagnosis.</p>	<p>Indications for Use of the modified device maintains the same general purpose and function as the predicate device, differing only by target MRI system.</p>

**510(k) SUMMARY**

<b>Characteristics</b>	<b>Predicate Device: InkSpace Imaging Body Array (K223487)</b>	<b>Subject Device: InkSpace Imaging Small Body Array</b>	<b>Substantial Equivalence Rationale</b>
	spine, shoulder, elbow, knee, foot, and prostate in GE 3.0T magnetic resonance imaging systems. These images, when interpreted by a trained physician, yield information that may assist in diagnosis.		
<b>Patient Anatomy and Population</b>	Imaging of general body part regions including cardiac regions; intended for pediatric patients. Imaging of general human anatomy, such as spine, shoulder, elbow, knee, foot, and prostate; intended for adult patients.	Imaging of general human anatomy in adult and pediatric patients, such as cardiac, spine, shoulder, elbow, knee, foot, and prostate.	N/A (Same).
<b>Comparison of Technological Characteristics</b>	24 channel, receive-only phased array, designed for 3.0T MR systems (127.73 MHz, hydrogen). Device is reusable, non-sterile, and by prescription only. Device is made of soft, pliable materials, and is composed of an anterior and a posterior pad.	24 channel, receive-only phased array, designed for 3.0T MR systems (123.2 MHz, hydrogen). Device is reusable, non-sterile, and by prescription only. Device is made of soft, pliable materials, and is composed of two identical 12ch pads.	N/A (Same).
<b>Compatible MRI Systems</b>	GE 3.0T MRI Systems	Siemens 3.0T MRI Systems	Differed target MRI system vendor. Same magnetic field strength.
<b>Biocompatibility</b>	ISO 10993 testing	ISO 10993 testing	N/A (Same).

## 510(k) SUMMARY

---

### SUBSTANTIAL EQUIVALENCE

The InkSpace Imaging Small Body Array is substantially equivalent to the predicate device, InkSpace Imaging Body Array with regard to function, intended use, and physical characteristics. There are no differences in the technological characteristics between the devices which would raise any different questions of safety or efficacy. Thus, the InkSpace Imaging Small Body Array is substantially equivalent to the predicate device.

### PERFORMANCE DATA [807.92(B)]

Performance testing was performed in accordance with the FDA Guidance Document entitled “*Magnetic Resonance (MR) Receive-only Coil – Performance Criteria for Safety and Performance Based Pathway*,” issued December 11, 2020.

Some of the nonclinical testing conducted on the InkSpace Imaging Body Array (predicate device) is still applicable to the InkSpace Imaging Small Body Array.

Additional analysis was conducted based on variations between the devices and differences between MR system vendors. These results and images listed below have been included as part of this submission.

### [807.92(b)(1)] Nonclinical Testing Summary:

The nonclinical, bench tests, per the FDA Guidance Document entitled “*Magnetic Resonance (MR) Receive-only Coil – Performance Criteria for Safety and Performance Based Pathway*,” were performed on the InkSpace Imaging Small Body Array:

- Image Signal to Noise (SNR)
- Image Uniformity
- Surface Heating
- Decoupling Circuit
- EMC – Immunity, Electrostatic Discharge
- General Electrical/Mechanical Safety
- Acquired Image Quality
- Biocompatibility Testing (Cytotoxicity, Sensitization, Irritation)

In addition to the above nonclinical testing, a risk analysis was performed to validate that the following testing performed on the InkSpace Imaging Body Array (predicate device) could be applied to the InkSpace Imaging Small Body Array without retesting:

- Biocompatibility Testing (Cytotoxicity, Sensitization, Irritation)
- Usability Testing
- Transit Testing

## **510(k) SUMMARY**

---

The collective results of the nonclinical testing demonstrate that the materials chosen, the manufacturing processes, and design of the InkSpace Imaging Small Body Array meet the established specifications necessary for consistent performance during its intended use. In addition, the collective bench testing demonstrates that the InkSpace Imaging Small Body Array does not raise different questions of safety or effectiveness for MR imaging examinations with Siemens 3.0T MRI Systems when compared to the predicate device.

### **[807.92(b)(2)] Clinical Testing Summary:**

Not applicable. Clinical testing was not required to demonstrate substantial equivalence of InkSpace Imaging Small Body Array to the predicate device. Acquired image quality was analyzed based on sample clinical images. These images, as well as a review of their clinical quality by a board-certified radiologist, have been included as part of this submission.

### **CONCLUSIONS [807.92(B)(3)]**

Based on the results from the nonclinical tests performed in support of the InkSpace Imaging Small Body Array, it is concluded that the modified device is safe, effective, and performs at least as safely and effectively as the legally marketed predicate device.

### **SUMMARY**

The InkSpace Imaging Small Body Array is substantially equivalent to the predicate device.