STOOL COLLECTION KIT FOR GI MOLECULAR TEST- container, specimen, non-sterile H&M MEDICAL COLLECTION KITS LLC

Stool collection kit for GI molecular test

KIT Includes

- Lab Tissue Mailer, Overall Height 3-7/8 In., Overall Length 4-1/8 In., Zoro Number: G3670484, Mfr Number: 729/CP; Qty-1
- "Mini Ice Pack from an Ice pack Sqaure Mat, Length 16-1/2 In., Width 14-3/4 In., Height 14-3/4 In., Zoro #: G5830133 | Mfr #: 405SQ (each matt has 60 pieces, which only one small piece is used)
- **Mailing Carton,** Capacity 72 cu. in., Color White, Inside Width 4 In., Inside Length 6 In., Inside Depth 3 In,Zoro#: G1072556 | Mfr#: 0101-1681010; Qty-1
- **One 3ml transfer pipette;** Transfer Pipet, 3.0mL, Small Bulb, Graduated to 1mL, 140mm, Bulb Draw 2.2mL, globe scientific, item#137035; Qty-1
- One pair of medium latex free gloves
- **Bag, Biohazard Specimen Transport**, 6" x 9", Ziplock with Document Pouch, Globe scientific,item#4919; Qty-1
- **Container, Fecal, 30mL,** attached Screw Cap with Spoon, PP, Conical Bottom, Self-Standing, globe scientific, item#109120; Qty-1
- **Sterile cotton tipped applicators**, Dynarexcorporation, 2 X 6"inch, reorder no. 4305, lot: 25519; Qty-1
- **Microtube,** 1.5mL, Self-Standing, Attached Screw Cap for Color Insert, with O-Ring, STERILE, PP, globe scientific, item#111720; Qty-2
- **Alcohol pads:** PT# 5110 PT# # 5110- Prep Surgical Alcohol Pads Webcol Large by, Kendall Company Qty: 2 in each kit
- **1A Mini zip lock with a toilet cover seat:** The Hospeco Health Gards half-fold toilet seat cover (part number:HG-1000) which is located inside a mini zip lock. Qty:1
- PURELL Sanitizing Hand Wipes Individually Wrapped: Qty;1 in each kit
- Labels for patient identification
- Instruction sheet for collection the sample
- Requisition form for the clinicians

Drug Facts

Active ingredient

Ethyl alcohol 62%

Purpose

Antimicrobial

Keep out of reach of children. If swallowed, get medical help or contact a Poison Control Center right away.

Use Hand sanitizer to help reduce bacteria on the skin

Warnings

- Flammable. Keep away from fire or flame,
- For external use only

- When using this product do not use in or near the eyes. In case of contact, rinse eyes thoroughly with water.
- **Stop use and ask a doctor if** irritation or rash appears and lasts

Directions • Unfold towelette; Wipe hands; Discard after use

Other information

- Do not store above 110°F (43°C)
- May discolor certain fabrics or surfaces

Inactive ingredients

Water (Aqua), Isopropyl Alcohol, Acetamidopropyl Trimonium Chloride, Aloe Barbadensis Leaf Juice, C12-15 Alkyl Benzoate, Cocamidopropyl PG-Dimonium Chloride Phosphate, Glycerin, Isopropyl Myristate, Propylene Glycol, Retinyl Palmitate, Tocopheryl Acetate, Zea Mays (Corn) Oil, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Aminomethyl Propanol, Fragrance (Parfum), Methylparaben

Questions?

1-800-321-9647

PACKAGE INFORMATION

NDC 21749-363-76

Purell Brand®

SANITIZING HAND WIPE

Kills 99.9% of most common germs that may make you sick

Not for retail sale

ONE WIPE

Distributed by: GOJO Industries, Inc. Akron OH 44309 Purell is a trademark of Warner-Lambert Company LLC

Packaging



Active ingredient

Active ingredient

Ethyl Alcohol 62%...

Antimicrobia

Use Hand sanitizer to help reduce bacteria on the skin

Warnings • Flammable. Keep away from fire or
flame. For external use only • When using this
product do not use in or near the eyes. In case of contact,
rinse eyes thoroughly with water. • Stop use and ask a
doctor if irritation or rash appears and lasts • Keep out
of reach of children. If swallowed, get medical help or
contact a Poison Control Center right away.

Directions Wipe hands, discard

Other information May discolor certain fabrics or surfaces

STOOL COLLECTION KIT FOR GI MOLECULAR TEST

container, specimen, non-sterile kit

Product Information

Product Type MEDICAL DEVICE Item Code (Source) NHRIC:70049-211

Packaging

l	# Item Code Package Description		Marketing Start Date	Marketing End Date	
ı	1 1	NHRIC:70049-211-70	1 in 1 KIT; Type 1: Convenience Kit of Co-Package		

Quantity of Parts

Part #	Package Quantity	Total Product Quantity
Part 1	1 APPLICATOR	1.75 mL

Part 1 of 1

PURELL SANITIZING HAND WIPES

ethyl alcohol cloth

T		C	. •
UMUU	luct Ir	ıforma	ation
1100	исти	HULLIN	1 LIVII

Item Code (Source)	NDC:21749-363
Route of Administration	TOPICAL

	Active Ingredient/Active Moiety				
ı	Ingredient Name	Basis of Strength	Strength		
ı	ALCOHOL (UNII: 3K9958V90M) (ALCOHOL - UNII:3K9958V90M)	ALCOHOL	0.62 mL in 1 mL		

Inactive Ingredients		
Ingredient Name	Strength	
WATER (UNII: 059QF0KO0R)		
ISOPROPYL ALCOHOL (UNII: ND2M416302)		
ALOE VERA LEAF (UNII: ZY8 1Z8 3H0 X)		
C12-15 ALKYL BENZO ATE (UNII: A9 EJ3J6 1HQ)		
GLYCERIN (UNII: PDC6A3C0OX)		
ISOPROPYL MYRISTATE (UNII: 0 RE8 K4LNJS)		
PROPYLENE GLYCOL (UNII: 6DC9Q167V3)		
VITAMIN A PALMITATE (UNII: 1D1K0 N0 VVC)		
ALPHA-TO CO PHERO L ACETATE (UNII: 9E8 X80 D2L0)		
CORN OIL (UNII: 8470G57WFM)		
AMINO METHYLPRO PANO L (UNII: LU49 E6626Q)		
METHYLPARABEN (UNII: A218 C7H19 T)		

Packaging					
# Item Code	Package Description	Marketing Start Date	Marketing End Date		
1 NDC:21749-363- 76	1 in 1 PACKET				
1	1.75 mL in 1 APPLICATOR; Type 0: Not a Combination Product				

Marketing Information					
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date		
OTC monograph not final	part333A	01/22/2007			

Marketing Information				
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date	
exempt device	NNI	09/01/2015		

Labeler - H&M MEDICAL COLLECTION KITS LLC (079931322)

Establishment			
Name	Address	ID/FEI	Business Operations
H&M MEDICAL COLLECTION KITS LLC		079931322	repack

Establishment					
Name	Address	ID/FEI	Business Operations		
GOJO Industries, Inc		036424534	manufacture		

Revised: 9/2015

H&M MEDICAL COLLECTION KITS LLC