

EASY GEL- jelly, lubricating, for transurethral surgical instrument gel
Leosons Overseas Corp

EASY GEL™

Use

For lubrication to provide easy insertion of catheters, endoscopes, or gloved fingers into bodily orifices.

Directions

Before first use, remove or puncture seal under cap. Squeeze lubricating jelly out of tube as needed. Discard after use.

Warnings

If irritation or discomfort occurs, discontinue use and consult a doctor. Keep out of reach of children. Avoid contact with eyes.

Other Information

Store at room temperature: 15–30° C (59–86° F). Lot number and expiration date on crimp. Sterility of contents guaranteed unless package is opened, damaged, or the seal is not intact.

Chlorhexidine Gluconate (CHG) was not used in the manufacture of this product.

PRINCIPAL DISPLAY PANEL - 113.4 G Tube Box

EASY•GEL™

STERILE

LUBRICATING JELLY

Water Soluble • Bacteriostatic

NET WT. 4 OZ (113.4 G)

Use: For lubrication to provide easy insertion of catheters, endoscopes, or gloved fingers into bodily orifices.

Directions: Before first use, remove or puncture seal under cap. Squeeze lubricating jelly out of tube as needed. Discard after use.

Warnings: If irritation or discomfort occurs, discontinue use and consult a doctor. Keep out of reach of children. Avoid contact with eyes.

Other Information: Store at room temperature: 15–30° C (59–86° F). Lot number and expiration date on crimp. Sterility of contents guaranteed unless package is opened, damaged, or the seal is not intact.

Chlorhexidine Gluconate (CHG) was not used in the manufacture of this product.

Manufactured for:
Leosons
10 Maryland Avenue,
Albany, New York 12205, USA
Call Toll Free +1-855-452-9500
info@Leosonsintl.com

MADE IN USA



إيزي جِل[™]

معقم

LUBRICATING JELLY

Water Soluble • Bacteriostatic

NET WT. 4 OZ (113.4 G)

EASY•GEL[™]

STERILE

LUBRICATING JELLY

Water Soluble • Bacteriostatic

NET WT. 4 OZ (113.4 G)


EASY•GEL[™]


STERILE

LUBRICATING JELLY

Water Soluble • Bacteriostatic

NET WT. 4 OZ (113.4 G)





EASY GEL

jelly, lubricating, for transurethral surgical instrument gel

Product Information

Product Type	MEDICAL DEVICE	Item Code (Source)	NHRIC:69626-0011
Route of Administration	RECTAL		

Other Ingredients

Ingredient Kind	Ingredient Name	Quantity
INGR	Water (UNII: 059QF0KO0R)	
INGR	Glycerin (UNII: PDC6A3C0OX)	
INGR	Carboxypolymethylene (UNII: 0A5MM307FC)	
INGR	Carboxymethylcellulose Sodium, Unspecified Form (UNII: K679OBS311)	
INGR	Methylparaben (UNII: A2I8C7HI9T)	
INGR	Edetate Disodium Anhydrous (UNII: 8NLQ36F6MM)	
INGR	Sodium Hydroxide (UNII: 55X04QC32I)	
INGR	Propylparaben (UNII: Z8IX2SC1OH)	

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:69626-0011-5	113.4 g in 1 TUBE; Type 0: Not a Combination Product		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
Exempt device	FHX	12/01/2016	

Labeler - Leosons Overseas Corp (148605470)**Establishment**

Name	Address	ID/FEI	Business Operations
Ultra Seal Corporation		085752004	MANUFACTURE

Revised: 12/2016

Leosons Overseas Corp