

GENERAL SURGERY KIT- general surgery tray (kit)
Centurion Medical Products

General Surgery Kit

DESCRIPTION

0.9% Sodium Chloride Injection, USP is sterile and nonpyrogenic. It is a parenteral solution containing sodium chloride in water for injection intended for intravenous administration. Each 100 ml of 0.9% sodium chloride injection, USP contains 900 mg sodium chloride in water for injection. The pH is 5.6 (4.5 to 7.0). This solution contains no bacteriostat, antimicrobial agent or added buffer and is intended only as a single-dose injection. When smaller doses are required, the unused portion should be discarded. 0.9% Sodium Chloride is a parenteral fluid and electrolyte replenisher and is a white crystalline powder freely soluble in water.

250 mL

NDC 0409-7983-25

VisIV™ Container

**0.9% SODIUM CHLORIDE
INJECTION, USP**

**EACH 100 mL CONTAINS SODIUM
CHLORIDE 900 mg IN WATER FOR
INJECTION. ELECTROLYTES PER 1000 mL:
SODIUM 154 mEq; CHLORIDE 154 mEq.
308 mOsmo/LITER (CALC).**

pH 5.6 (4.5 to 7.0)

**ADDITIVES MAY BE INCOMPATIBLE.
CONSULT WITH PHARMACIST, IF
AVAILABLE. WHEN INTRODUCING
ADDITIVES, USE ASEPTIC TECHNIQUE,
MIX THOROUGHLY AND DO NOT
STORE. SINGLE-DOSE CONTAINER. FOR
INTRAVENOUS USE. USUAL DOSAGE:
SEE INSERT. STERILE, NONPYROGENIC.
STORE AT 20 TO 25°C (68 TO 77°F).
IFEE USD CONTROLLED ROOM**

50

100

[SEE USP CONTROLLED ROOM TEMPERATURE.] PROTECT FROM FREEZING. SEE INSERT. USE ONLY IF SOLUTION IS CLEAR AND CONTAINER IS UNDAMAGED. MUST NOT BE USED IN SERIES CONNECTIONS. PVC-FREE, DEHP-FREE. VisiV IS A TRADEMARK OF HOSPIRA. DO NOT REMOVE CAPS UNTIL READY FOR USE. IF LEAKS ARE FOUND, DISCARD SOLUTION AS STERILITY MAY BE IMPAIRED.

150

200

Rx ONLY



IM-1693 (7/08)

PRINTED IN USA
HOSPIRA, INC.
LAKE FOREST, IL 60045 USA



MM

CENTURION®

032312G

Example General

Reorder

Surgery Tray

XXXXXXXXXX

CONTENTS:

1 STERILE EXAMPLE GENERAL SURGERY TRAY


1 STERILE SODIUM CHLORIDE SOLUTION 500mL

Example Only - Components & Title May Vary

LOT 0000000000



N/A



(01) 1 0653160 00000 6

NOTES:

Store between 20-25°C (68-77°F).



Do not reuse.
Single use only.

LATEX FREE



Consult
instructions
for use.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

MM

GENERAL SURGERY KIT

general surgery tray (kit) kit

Product Information

Product Type	MEDICAL DEVICE	Item Code (Source)	NHRIC:24840-1535
---------------------	----------------	---------------------------	------------------

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:24840-1535-2	10 in 1 CASE		
1	NHRIC:24840-1535-1	1 in 1 PACKAGE, COMBINATION		

Quantity of Parts

Part #	Package Quantity	Total Product Quantity
Part 1	1 BAG	500 mL

Part 1 of 1**SODIUM CHLORIDE**

sodium chloride injection, solution

Product Information

Item Code (Source)	NDC:0409-7983
---------------------------	---------------

Route of Administration	INTRAVENOUS
--------------------------------	-------------

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
SODIUM CHLORIDE (UNII: 451W47IQ8X) (SODIUM CATION - UNII:L4M0NH37, CHLORIDE ION - UNII:Q32ZN48698)	SODIUM CHLORIDE	900 mg in 100 mL

Inactive Ingredients

Ingredient Name	Strength
WATER (UNII: 059QF0KO0R)	

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
---	-----------	---------------------	----------------------	--------------------

1	NDC:0409-7983-30	500 mL in 1 BAG		
Marketing Information				
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date	
ANDA	ANDA016366	04/10/2006		
Marketing Information				
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date	
exempt device	LRO	01/01/2012		

Labeler - Centurion Medical Products (017246562)

Establishment			
Name	Address	ID/FEI	Business Operations
Centurion Medical Products		017246562	manufacture, repack

Establishment			
Name	Address	ID/FEI	Business Operations
Centurion Medical Products		148522279	manufacture, repack

Establishment			
Name	Address	ID/FEI	Business Operations
Centurion Medical Products		626660810	manufacture, repack

Establishment			
Name	Address	ID/FEI	Business Operations
Hospira Inc		093132819	manufacture