

CYGZYME PLUS- peracetic acid liquid Valley Inc.

Disclaimer: This drug has not been found by FDA to be safe and effective, and this labeling has not been approved by FDA. For further information about unapproved drugs, click here.

ACTIVE INGREDIENT

Peracetic acid 0.05%

INACTIVE INGREDIENTS

Water, Lauramine oxide, Sodium bicarbonate, Sodium acetate, anhydrous, Subtilisins, Cellulase, Rizolipase, alpha Amylase, Brilliant Blue, Alcohols, C12-14, ethoxylated, propoxylated, Sodium sesquicarbonate

PURPOSE

Antifungal

WARNINGS

- For external use only
- Do not use on the human body
- Avoid contact with eyes.
- Discontinue use if signs of irritation or rashes appear.
- Keep out of reach of children.

KEEP OUT OF REACH OF CHILDREN

Keep out of reach of children

Uses

Helps prevent sterilant

Directions

- Rinse thoroughly for more than 5 minutes. Rinse for more than 15 seconds with sterile purified water flowing out after soaking.
- After the disinfection of the medical instrument is completed, confirm the persistence.
- Excess water may be used depending on the purpose of use.

QUESTIONS

+82 02 1566 2723

PACKAGE LABEL.PRINCIPAL DISPLAY PANEL

<p>Must keep stand during storage & transportation</p>																													
<p>Cygzyme plus</p> <table border="1"> <tr> <th colspan="2">Drug Facts</th> </tr> <tr> <th>Active ingredient(s)</th> <th>Purpose</th> </tr> <tr> <td>Peracetic acid 0.05%</td> <td>Antifungal</td> </tr> <tr> <th colspan="2">Use(s)</th> </tr> <tr> <td colspan="2"> <ul style="list-style-type: none"> Helps prevent sterilant. </td> </tr> <tr> <th colspan="2">Warnings</th> </tr> <tr> <td colspan="2"> <ul style="list-style-type: none"> For external use only. Do not use on the human body. Avoid contact with eyes. Discontinue use if signs of irritation or rashes appear. Keep out of reach of children. </td> </tr> <tr> <th colspan="2">Directions</th> </tr> <tr> <td colspan="2"> <ul style="list-style-type: none"> Rinse thoroughly for more than 5 minutes. Rinse for more than 15 seconds with sterile purified water flowing out after soaking. After the disinfection of the medical instrument is completed, confirm the persistence. Excess water may be used depending on the purpose of use. </td> </tr> <tr> <th colspan="2">Questions</th> </tr> <tr> <td colspan="2"> <ul style="list-style-type: none"> +82 053 1566 2723 </td> </tr> <tr> <th colspan="2">Inactive ingredients</th> </tr> <tr> <td colspan="2"> Water, Lauramine oxide, Sodium bicarbonate, Sodium acetate, anhydrous, Stibilsins, Cellulase, Ricolipase, alpha Amylase, Brilliant Blue, Alcohols, C12-14 ethoxylated, propoxylated, Sodium sesquicarbonate </td> </tr> </table> <p>Distributor & Manufacturer: Valley Inc 16, Padong-ro 32-gil, Suseong-gu, Daegu, Republic of Korea 4,000ml / 135.26 Fl. Oz</p>	Drug Facts		Active ingredient(s)	Purpose	Peracetic acid 0.05%	Antifungal	Use(s)		<ul style="list-style-type: none"> Helps prevent sterilant. 		Warnings		<ul style="list-style-type: none"> For external use only. Do not use on the human body. Avoid contact with eyes. Discontinue use if signs of irritation or rashes appear. Keep out of reach of children. 		Directions		<ul style="list-style-type: none"> Rinse thoroughly for more than 5 minutes. Rinse for more than 15 seconds with sterile purified water flowing out after soaking. After the disinfection of the medical instrument is completed, confirm the persistence. Excess water may be used depending on the purpose of use. 		Questions		<ul style="list-style-type: none"> +82 053 1566 2723 		Inactive ingredients		Water, Lauramine oxide, Sodium bicarbonate, Sodium acetate, anhydrous, Stibilsins, Cellulase, Ricolipase, alpha Amylase, Brilliant Blue, Alcohols, C12-14 ethoxylated, propoxylated, Sodium sesquicarbonate		 <p>HANDLE WITH CARE</p>	<p>Cygzyme plus</p> <p>4,000ml / 135.26 Fl. Oz</p>	<p>Cleaner and disinfectant for endoscopes & other medical instruments</p>
Drug Facts																													
Active ingredient(s)	Purpose																												
Peracetic acid 0.05%	Antifungal																												
Use(s)																													
<ul style="list-style-type: none"> Helps prevent sterilant. 																													
Warnings																													
<ul style="list-style-type: none"> For external use only. Do not use on the human body. Avoid contact with eyes. Discontinue use if signs of irritation or rashes appear. Keep out of reach of children. 																													
Directions																													
<ul style="list-style-type: none"> Rinse thoroughly for more than 5 minutes. Rinse for more than 15 seconds with sterile purified water flowing out after soaking. After the disinfection of the medical instrument is completed, confirm the persistence. Excess water may be used depending on the purpose of use. 																													
Questions																													
<ul style="list-style-type: none"> +82 053 1566 2723 																													
Inactive ingredients																													
Water, Lauramine oxide, Sodium bicarbonate, Sodium acetate, anhydrous, Stibilsins, Cellulase, Ricolipase, alpha Amylase, Brilliant Blue, Alcohols, C12-14 ethoxylated, propoxylated, Sodium sesquicarbonate																													
<p>Manufacturer VALLEY INC</p>																													

CYGZYME PLUS

peracetic acid liquid

Product Information

Product Type	HUMAN OTC DRUG	Item Code (Source)	NDC:81803-050
Route of Administration	TOPICAL		

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
PERACETIC ACID (UNII: I6KPI2E1HD) (PERACETIC ACID - UNII:I6KPI2E1HD)	PERACETIC ACID	0.05 g in 100 mL

Inactive Ingredients

Ingredient Name	Strength
WATER (UNII: 059QF0KO0R)	
LAURAMINE OXIDE (UNII: 4F6FC4MI8W)	

SODIUM BICARBONATE (UNII: 8MDF5V39QO)

SODIUM ACETATE (UNII: 4550K0SC9B)

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:81803-050-02	1 in 1 CARTON	02/01/2024	
1	NDC:81803-050-01	4000 mL in 1 BOTTLE; Type 0: Not a Combination Product		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
unapproved drug other		02/01/2024	

Labeler - Valley Inc. (695095445)

Registrant - Valley Inc. (695095445)

Establishment

Name	Address	ID/FEI	Business Operations
Valley Inc.		695095445	manufacture(81803-050)

Revised: 12/2024

Valley Inc.