

**ATLAS PREGNANCY TEST- visual, pregnancy hcg, prescription use
ECLIPSE PHARMA LLC**

Atlas Pregnancy Test

WARNINGS

- Check if the pouch is completely sealed.
- Keep away from children.
- Store at 36-86°F
- Preg Product



Product Packaging

Product Packaging




Preg cont back


Atlas Home Pregnancy Test

Use Instructions

- Remove the cap from the device, and place the absorbent tip directly into the urine stream for at least 10 seconds.
- Read results after 3 minutes.



Results



- Keep away from children.
- Store at 36-86°F.

Pack Content:

- Test Device.
- Instructions for use.

- Atlas Pregnancy Test
- Atlas Pregnancy Test
- customer.service@AtlasHomeTests.com
- www.AtlasHomeTests.com/en|Pregnancy_m

Atlas Medical Product Helpline
 Monday - Friday 9:00a.m. - 5:00 p.m.
 1-888-257-6332
 customer.service@atlashometests.com

Atlas Medical Product Helpline
 Monday - Friday 9:00 a.m. - 5:00 p.m.
 1-888-257-6332
 customer.service@atlashometests.com

Manufactured for Atlas Medical GmbH
 Blankenfelde-Mahlow, Germany
 www.atlas-medical.com

Atlas Medical

CE 0459

REF 70172001

ST IVD

36°F - 86°F



Step by step
Instructions video

ATLAS PREGNANCY TEST

visual, pregnancy hcg, prescription use

Product Information

Product Type	PRESCRIPTION MEDICAL DEVICE	Item Code (Source)	NHRIC:77468-200
---------------------	-----------------------------	---------------------------	-----------------

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:77468-200-01	1 in 1 PACKET		
1		1 in 1 KIT; Type 1: Convenience Kit of Co-Package		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
premarket notification	K102760	11/10/2023	

Labeler - ECLIPSE PHARMA LLC (108358472)

Revised: 11/2023

ECLIPSE PHARMA LLC