

LIDOBAND- lidocaine ring
Alberta Veterinary Laboratories Ltd

Disclaimer: This drug has not been found by FDA to be safe and effective, and this labeling has not been approved by FDA. For further information about unapproved drugs, click here.

LidoBand Castration Bands

LidoBand Castration Bands

Solvat

LidoBand Castration Bands

Each band contains 80 mg of lidocaine.

For castration of calves and lambs, and to deliver lidocaine as a local anesthetic. For use in calves under 250 pounds or lambs under 50 pounds.

FOR ANIMAL USE ONLY

KEEP OUT OF REACH OF CHILDREN

TAKE TIME OBSERVE LABEL DIRECTIONS

Description:

LidoBand™ is a latex rubber band designed to be placed around the base of the scrotum of calves and lambs for castration. Each castration band is infused with 80 mg of lidocaine, which provides local anesthesia to the application area beginning at 2 hours post-application and lasting up to 42 days.

Precautions:

Lidocaine is typically well tolerated. As with all local anesthetics, untoward effects may occur due to hypersensitivity.

LidoBand should not be used on animals with a known lidocaine hypersensitivity.

User Safety:

If skin comes in contact with LidoBand, wash affected area immediately with soap and water.

Storage:

Keep stored in a cool, dry place.

Do not re-use the band.

Directions for use: (8, 40, or 80-count box)

Each band is individually packaged in a foil-sealed blister. Do not touch the band with bare fingers. Wear gloves when applying LidoBand. Only apply LidoBand with a commercial band applicator (not included).

Step 1. Puncture the blister surface with the band applicator.

Step 2. Remove the band from packaging while loaded on the prongs of the applicator.

Step 3. Squeeze the applicator handle to spread the band.

Step 4. Pull the testicles into the scrotum, place the band above the testicles, and remove the applicator.

Directions for use: (250 count jar)

Bands contained within are to be used with a hopper chamber (not included). Do not touch the band with bare fingers. Wear gloves when applying LidoBand. Only apply LidoBand with a commercial band applicator (not included). Pour unused bands back into jar.

Step 1. Put bands into hopper chamber, and rotate to direct band into dispenser.

Step 2. Stick the prongs of the applicator into the band loop and remove from dispenser.

Step 3. Squeeze the applicator handle to spread the band.

Step 4. Pull the testicles into the scrotum, place the band above the testicles, and remove the applicator.

Product Labels



lidoband™ solvet

Castration Bands

Each band contains 80 mg of lidocaine.

For castration of calves and lambs, and to deliver lidocaine as a local anesthetic. For use in calves under 250 pounds or lambs under 50 pounds.

Foil
LidoBand
Blister

TAKE TIME
 OBSERVE LABEL DIRECTIONS

FOR ANIMAL USE ONLY
KEEP OUT OF REACH OF CHILDREN
NDC: 86174-101-03

40 Castration Bands

lidoband™ solvet

Castration Bands

Step 1. Puncture the blister surface with the band applicator.

Step 2. Remove the band from packaging while loaded on the prongs of the applicator.

Step 3. Squeeze the applicator handle to spread the band.

Step 4. Pull the testicles into the scrotum, place the band above the testicles, and remove the applicator.

Manufactured by :
AVL/Solvet Ltd. 7226 107 Ave, Calgary, AB T2C 6N6
(877) 456-2755. www.solvet.ca

lidoband™ solvet

Castration Bands

For castration of calves and lambs, and to deliver lidocaine as a local anesthetic. For use in calves under 250 pounds or lambs under 50 pounds.

Description: LidoBand™ is a latex rubber band designed to be placed around the base of the scrotum of calves and lambs for castration. Each castration band is infused with 80 mg of lidocaine, which provides local anesthesia beginning at 2 hours post-application and lasting up to 42 days.

Directions for use:
Each band is individually packaged in a foil-sealed blister. Do not touch the band with bare fingers. Wear gloves when applying LidoBand. Only apply LidoBand with a commercial band applicator (not included).

California requires:
WARNING: Cancer Risk
www.P65Warnings.ca.gov

U.S. Patent #: 11,596,510
Solvet is a trademark of Alberta Veterinary Laboratories Ltd.

lidoband™ solvet

Castration Bands

Precautions: Lidocaine is typically well tolerated. As with all local anesthetics, untoward effects may occur due to hypersensitivity. LidoBand should not be used on animals with a known lidocaine hypersensitivity.

User Safety: If skin comes in contact with LidoBand, wash affected area immediately with soap and water.

Storage: Keep stored in a cool, dry place. Do not re-use the band.

Lot:
Exp:

Manufactured by :
AVL/Solvet Ltd. 7226 107 Ave, Calgary, AB T2C 6N6
(877) 456-2755. www.solvet.ca

LidoBand 40-Count Box (8 bands per blister pack, 5 blister packs per box)

80 Castration Bands

FOR ANIMAL USE ONLY
KEEP OUT OF REACH OF CHILDREN

OBSERVE LABEL DIRECTIONS

TAKE TIME

Foil
LidoBand
Blister

Each band contains 80 mg of lidocaine.
Castration Bands
LidoBand™
NDC 86174-101-01

For castration of calves and lambs, and to deliver lidocaine as a local anesthetic. For use in calves under 250 pounds or lambs under 50 pounds.

lidoband solvet

Castration Bands
Each band contains 80 mg of lidocaine.

For castration of calves and lambs, and to deliver lidocaine as a local anesthetic. For use in calves under 250 pounds or lambs under 50 pounds.

Foil
LidoBand
Blister

FOR ANIMAL USE ONLY
KEEP OUT OF REACH OF CHILDREN
NDC: 86174-101-01

TAKE TIME
OBSERVE LABEL DIRECTIONS

80 Castration Bands

lidoband

Castration Bands

Directions for use:
Each band is individually packaged in a foil-sealed blister. Do not touch the band with bare fingers. Wear gloves when applying LidoBand. Only apply LidoBand with a commercial band applicator (not included).

Step 1. Push the blister surface with the band applicator.
Step 2. Remove the band from packaging while loaded on the prongs of the applicator.
Step 3. Separate the applicator handle to spread the band.
Step 4. Pull the testicles into the scrotum, place the band above the testicles, and remove the applicator.

80 Castration Bands

Manufactured by :
AVL/Solvet Ltd. 7226 107 Ave, Calgary, AB T2C 5N6
(877) 456-2755. www.solvet.ca

lidoband solvet

Castration Bands

For castration of calves and lambs, and to deliver lidocaine as a local anesthetic. For use in calves under 250 pounds or lambs under 50 pounds.

Description: LidoBand™ is a latex rubber band designed to be placed around the base of the scrotum of calves and lambs for castration. Each castration band is infused with 80 mg of lidocaine, which provides local anesthesia beginning at 2 hours post-application and lasting up to 42 days.

Precautions: Lidocaine is typically well tolerated. As with all local anesthetics, untoward effects may occur due to hypersensitivity. LidoBand should not be used on animals with a known lidocaine hypersensitivity.

User Safety: If skin comes in contact with LidoBand, wash affected area immediately with soap and water.

Storage: Keep stored in a cool, dry place. Do not re-use the band.

California requires:
WARNING: Cancer Risk
www.P65Warnings.ca.gov

80 Castration Bands

U.S. Patent #: 11,596,510
Solvet is a trademark of Alberta Veterinary Laboratories Ltd.

lidoband

Castration Bands

Directions for use:
Each band is individually packaged in a foil-sealed blister. Do not touch the band with bare fingers. Wear gloves when applying LidoBand. Only apply LidoBand with a commercial band applicator (not included).

Step 1. Push the blister surface with the band applicator.
Step 2. Remove the band from packaging while loaded on the prongs of the applicator.
Step 3. Separate the applicator handle to spread the band.
Step 4. Pull the testicles into the scrotum, place the band above the testicles, and remove the applicator.

80 Castration Bands

Manufactured by :
AVL/Solvet Ltd. 7226 107 Ave, Calgary, AB T2C 5N6
(877) 456-2755. www.solvet.ca

LidoBand 80-Count Box (8 bands per blister pack, 10 blister packs per box)

Description: LidoBand™ is a latex rubber band designed to be placed around the base of the scrotum of calves and lambs for castration. Each castration band is infused with 80 mg of lidocaine, which provides local anesthesia to the application area beginning at 2 hours post-application and lasting up to 42 days.

Precautions: Lidocaine is typically well tolerated. As with all local anesthetics, untoward effects may occur due to hypersensitivity. LidoBand should not be used on animals with a known lidocaine hypersensitivity.

User Safety: If skin comes in contact with LidoBand, wash affected area immediately with soap and water.

Storage: Keep stored in a cool, dry place. Do not re-use the band.

TAKE TIME
OBSERVE LABEL DIRECTIONS

U.S. Patent #: 11,596,510
Manufactured by :
AVL/Solvet Ltd. 7226 107 Ave, Calgary, AB T2C 5N6
(877) 456-2755. www.solvet.ca
Solvet is a trademark of Alberta Veterinary Laboratories Ltd.

solvet NDC: 86174-101-02

lidoband

Castration Bands
Each band contains 80 mg of lidocaine

For castration of calves and lambs, and to deliver lidocaine as a local anesthetic. For use in calves under 250 pounds or lambs under 50 pounds.

FOR ANIMAL USE ONLY
KEEP OUT OF REACH OF CHILDREN

250 castration bands

Directions for use:
Bands contained within are to be used with a hopper chamber (not included). Do not touch the band with bare fingers. Wear gloves when applying LidoBand. Only apply LidoBand with a commercial band applicator (not included). Pour unused bands back into jar.

Step 1. Put bands into hopper chamber, and rotate to direct band into dispenser.
Step 2. Stick the prongs of the applicator into the band loop and remove from dispenser.
Step 3. Squeeze the applicator handle to spread the band.
Step 4. Pull the testicles into the scrotum, place the band above the testicles, and remove the applicator.

LOT :
EXP. :

California requires:
WARNING: Cancer Risk
www.P65Warnings.ca.gov

LidoBand 250-Count Jar

SEE INSIDE COVER

DIRECTIONS FOR USE

NOT FOR RETAIL SALE

www.lidoband.com

U.S. Patent #: 11,596,510

lidoband™ solvet

Castration Bands


Each band contains 80 mg of lidocaine.

For castration of calves and lambs, and to deliver lidocaine as a local anesthetic. For use in calves under 250 pounds or lambs under 50 pounds.

FOR ANIMAL USE ONLY

KEEP OUT OF REACH OF CHILDREN

NDC: 86174-101-04



8 Castration Bands

www.lidoband.com

Manufactured by :
AVL/Solvet Ltd. 7226 107 Ave, Calgary, AB T2C 5N6
(877) 456-2755. www.solvet.ca

Description: LidoBand™ is a latex rubber band designed to be placed around the base of the scrotum of calves and lambs for castration. Each castration band is infused with 80 mg of lidocaine, which provides local anesthesia beginning at 2 hours post-application and lasting up to 42 days.

Directions for use:

Each band is individually packaged in a foil-sealed blister. Do not touch the band with bare fingers. Wear gloves when applying LidoBand. Only apply LidoBand with a commercial band applicator (not included).

Precautions: Lidocaine is typically well tolerated. As with all local anesthetics, untoward effects may occur due to hypersensitivity. LidoBand should not be used on animals with a known lidocaine hypersensitivity.

User Safety: If skin comes in contact with LidoBand, wash affected area immediately with soap and water.


Storage: Keep stored in a cool, dry place.
Do not re-use the band.

www.lidoband.com




California requires:
WARNING: Cancer Risk
www.P65Warnings.ca.gov


lidoband[™] solvet




Step 1. Puncture the blister surface with the band applicator.



Step 2. Remove the band from packaging while loaded on the prongs of the applicator.



Step 3. Squeeze the applicator handle to spread the band.



Step 4. Pull the testicles into the scrotum, place the band above the testicles, and remove the applicator.

LidoBand 8-Count Box (8 bands per blister pack, 1 blister pack per box)

LIDOBAND

lidocaine ring

Product Information

Product Type	OTC ANIMAL DRUG	Item Code (Source)	NDC:86174-101
Route of Administration	TOPICAL		

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
LIDOCAINE (UNII: 98PI200987) (LIDOCAINE - UNII:98PI200987)	LIDOCAINE	80 mg

Inactive Ingredients

Ingredient Name	Strength
ISOPROPYL MYRISTATE (UNII: 0RE8K4LNJS)	30 mg

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:86174-101-01	80 in 1 BOX		
2	NDC:86174-101-02	250 in 1 JAR		
3	NDC:86174-101-03	40 in 1 BOX		
4	NDC:86174-101-04	8 in 1 BOX		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
unapproved drug other		06/20/2023	

Labeler - Alberta Veterinary Laboratories Ltd (204293518)

Registrant - Alberta Veterinary Laboratories Ltd (204293518)

Establishment

Name	Address	ID/FEI	Business Operations
Alberta Veterinary Laboratories Ltd		204293518	manufacture

Establishment

Name	Address	ID/FEI	Business Operations
Swati Spentose Private Limited		650621472	api manufacture

Revised: 9/2023

Alberta Veterinary Laboratories Ltd