

X-STATIC SILVERSEAL HYDROGEL DRESSING 4 X 5- dressing, wound and burn, hydrogel w/drug and/or biologic
X-STATIC SILVERSEAL HYDROGEL DRESSING 2 X 3- dressing, wound and burn, hydrogel w/drug and/or biologic
X-STATIC SILVERSEAL HYDROGEL DRESSING 2 X 2- dressing, wound and burn, hydrogel w/drug and/or biologic
X-STATIC SILVERSEAL HYDROGEL DRESSING 3- dressing, wound and burn, hydrogel w/drug and/or biologic
HYDRESS- dressing,wound,hydrophilic
Alliqua Inc

HYDROGEL DRESSING

Apply as a single layer.

Overlay with conventional dressings.

SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier

Silverseal® Hydrogel is intended for the management of wounds and to provide an antimicrobial barrier

Applicable to:

Partial and full thickness ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.

Silverseal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber

Antibacterial activity was confirmed via in vitro studies of bacterial growth dressing.

To facilitate removal, moisten with sterile water.

Caution: federal law restricts this device to sale by or on the order at a licensed practitioner.

Use before expiration date.

Principal Display Panel – 2” X 2” Dressing Box Label

SILVER
SEAL®

(00)00000000000000

HYDROGEL DRESSING

NDC 58198-122-10 PRODUCT CODE

AA-000000-000

STERILE Contents: Ten (10) Dressing
SINGLE USE ONLY 2” X 2” (5.1cm X 5.1cm)

Manufactured by Aquamed Technologies Inc.
For Alliqua Biomedical Inc.

2150 Cabot Boulevard West, Langhorne, PA 19047
www.Alliqua.com



**SILVER
SEAL®**



(00)00000000000000

HYDROGEL DRESSING

PRODUCT CODE

AA-000000-000

NDC 58198-122-10

**STERILE
SINGLE USE ONLY**

**Contents: Ten (10) Dressing
2" X 2" (5.1cm X 5.1cm)**



Apply as a single layer.
Overlay with conventional dressings.
SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier



Silverseal® Hydrogel is intended for the management of wounds and to provide an antimicrobial barrier



Applicable to:
Partial and full thickness ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.
Silverseal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber

Antibacterial activity was confirmed via in vitro studies of bacterial growth dressing.

To facilitate removal, moisten with sterile water.

Caution: federal law restricts this device to sale by or on the order at a licensed practitioner.

Use before expiration date.

Manufactured by Aquamed Technologies Inc.
For Alliqua Biomedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047
www.Alliqua.com

Principal Display Panel – 2” X 2” Dressing Pouch Label

SILVER

SEAL[®]

(00)00000000000000

HYDROGEL DRESSING

NDC 58198-122-01 PRODUCT CODE

AA-000000-000

STERILE Contents: One (1) Dressing

SINGLE USE ONLY 2” X 2” (5.1cm X 5.1cm)

Manufactured by Aquamed Technologies Inc.

For Alliqua Biomedical Inc.

2150 Cabot Boulevard West, Langhorne, PA 19047

www.Alliqua.com



**SILVER
SEAL®**



(00)00000000000000

HYDROGEL DRESSING

NDC 58198-122-01

**STERILE
SINGLE USE ONLY**

PRODUCT CODE

AA-000000-000

**Contents: One (1) Dressing
2" X 2" (5.1cm X 5.1cm)**



Apply as a single layer.
Overlay with conventional dressings.
SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier



Silverseal® Hydrogel is intended for the management of wounds and to provide an antimicrobial barrier



Applicable to:
Partial and full thickness ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.
Silverseal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber

Antibacterial activity was confirmed via in vitro studies of bacterial growth dressing.

To facilitate removal, moisten with sterile water.

Caution: federal law restricts this device to sale by or on the order at a licensed practitioner.

Use before expiration date.

Manufactured by Aquamed Technologies Inc.
For Alliqua Biomedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047

Principal Display Panel – 2” X 3” Dressing Box Label

NDC 58198-123-10

SILVER
SEAL[®]

(01)10816789010027

HYDROGEL DRESSING

PRODUCT CODE
AB-200300-103

STERILE Contents: Ten (10) dressings
SINGLE USE ONLY 2” X 3” (5.1cm X 7.6cm)

Manufactured by Aquamed Technologies Inc.


for Alliqua Biomedical Inc.

2150 Cabot Boulevard West, Langhorne, PA 19047

www.Alliqua.com

**SILVER
SEAL®**

NDC#58189-123-10



(01)10816789010027

HYDROGEL DRESSING

PRODUCT CODE

AB-200300-103



STERILE

SINGLE USE ONLY

Contents: Ten (10) dressings

2" X 3" (5.1cm X 7.6cm)



Apply as a single layer. Overlay with conventional dressings. SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier. To facilitate removal, moisten with sterile water.



Applicable to: partial and full thickness dermal ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.



SilverSeal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber. Antibacterial activity was confirmed via in vitro studies of bacterial growth in the dressing.

Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.
Use before expiration date.

Manufactured by AquaMed Technologies Inc.
for Alliqua BioMedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047
www.Alliqua.com

Principal Display Panel – 4" X 5" Dressing Box Label

NDC 58198-123-01

SILVER
SEAL®

(01)00816789010020 HYDROGEL DRESSING

PRODUCT CODE

AB-200300-103

STERILE Contents: One (1) dressing
SINGLE USE ONLY 2" X 3" (5.1cm X 7.6cm)

Manufactured by Aquamed Technologies Inc.

for Alliqua Biomedical Inc.

2150 Cabot Boulevard West, Langhorne, PA 19047

NDC#58189-123-01



(01) 00816789010020

HYDROGEL DRESSING

PRODUCT CODE
AB-200300-103

.....
STERILE
SINGLE USE ONLY

Contents: One (1) dressing
2" X 3" (5.1cm X 7.6cm)



Apply as a single layer. Overlay with conventional dressings. SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier. To facilitate removal, moisten with sterile water.

Applicable to: partial and full thickness dermal ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.

SilverSeal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber. Antibacterial activity was confirmed via in vitro studies of bacterial growth in the dressing.

Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner. Use before expiration date.

Manufactured by Aquamed Technologies Inc.
for Alliqua BioMedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047

Principal Display Panel – 3” Circle Dressing Box Label

SILVER

SEAL[®]

(00)00000000000000

HYDROGEL DRESSING

NDC 58198-130-10 PRODUCT CODE

AA-000000-000

STERILE Contents: Ten (10) Dressing
SINGLE USE ONLY 3” Circle (7.6cm circle)

Manufactured by Aquamed Technologies Inc.

For Alliqua Biomedical Inc.

2150 Cabot Boulevard West, Langhorne, PA 19047

www.Alliqua.com



**SILVER
SEAL®**



(00)00000000000000

HYDROGEL DRESSING

NDC 58198-130-10

**STERILE
SINGLE USE ONLY**

PRODUCT CODE

AA-000000-000

**Contents: Ten (10) Dressing
3" Circle (7.6cm circle)**



Apply as a single layer.
Overlay with conventional dressings.
SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier



Silverseal® Hydrogel is intended for the management of wounds and to provide an antimicrobial barrier



Applicable to:
Partial and full thickness ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.
Silverseal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber

Antibacterial activity was confirmed via in vitro studies of bacterial growth dressing.

To facilitate removal, moisten with sterile water.

Caution: federal law restricts this device to sale by or on the order at a licensed practitioner.

Use before expiration date.

Manufactured by Aquamed Technologies Inc.
For Alliqua Biomedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047
www.Alliqua.com

Principal Display Panel – 3” Circle Dressing Pouch Label

SILVER

SEAL[®]

(00)00000000000000

HYDROGEL DRESSING

NDC 58198-130-01 PRODUCT CODE

AA-000000-000

STERILE Contents: One (1) Dressing

SINGLE USE ONLY 3” Circle (7.6cm circle)

Manufactured by Aquamed Technologies Inc.

For Alliqua Biomedical Inc.

2150 Cabot Boulevard West, Langhorne, PA 19047

www.Alliqua.com



**SILVER
SEAL®**



(00)00000000000000

HYDROGEL DRESSING

PRODUCT CODE

AA-000000-000

NDC 58198-130-01

**STERILE
SINGLE USE ONLY**

**Contents: One (1) Dressing
3" Circle (7.6cm circle)**



Apply as a single layer.
Overlay with conventional dressings.
SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier



Silverseal® Hydrogel is intended for the management of wounds and to provide an antimicrobial barrier



Applicable to:
Partial and full thickness ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.
Silverseal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber

Antibacterial activity was confirmed via in vitro studies of bacterial growth dressing.

To facilitate removal, moisten with sterile water.

Caution: federal law restricts this device to sale by or on the order at a licensed practitioner.

Use before expiration date.

Manufactured by Aquamed Technologies Inc.
For Alliqua Biomedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047
www.Alliqua.com

Principal Display Panel – 4" X 4" Circle Dressing Box Label

NDC 58198-044-10 PRODUCT CODE

AB-400400-001

HYDRESSTM

(01)10816789010003

HYDROGEL DRESSING


STERILE Contents: Ten (10) Hydrogel Dressings

SINGLE USE ONLY 4" X 4" (10cm X 10cm)


Manufactured by Aquamed Technologies Inc. for Alliqua Biomedical Inc.

2150 Cabot Boulevard West, Langhorne, PA 19047 | Tel 855-425-5478 | www.Alliqua.com

NDC#58189-044-10



PRODUCT CODE
AB-400400-001



(01)10816789010003

HYDROGEL DRESSING

STERILE
SINGLE USE ONLY

Contents: Ten (10) Hydrogel Dressings
4" X 4" (10cm X 10cm)

Hydress™ should be placed on the wound and bandaged as necessary.



Directions: Remove dressing from packaging and remove the blue plastic sheet from dressing. Apply exposed side of dressing to wound and bandage as necessary. Change dressing daily.

For use on minor burns, minor cuts and scrapes.

Ingredients: Deionized water, polyethylene oxide, methylparaben, propylparaben and citric acid.

To facilitate removal, moisten with sterile water.

Use before expiration date.

Manufactured by AquaMed Technologies Inc. for Alliqua BioMedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047 | Tel 855-425-5478 | www.Alliqua.com

Principal Display Panel – 4" X 4" Circle Dressing Pouch Label

NDC 58198-044-10 PRODUCT CODE

AB-400400-001

HYDRESS_{TM}

(01)00816789010006

HYDROGEL DRESSING

STERILE Contents: One (1) dressing

SINGLE USE ONLY 4" X 4" (10cm X 10cm)


Manufactured by Aquamed Technologies Inc. for Alliqua Biomedical Inc.

2150 Cabot Boulevard West, Langhorne, PA 19047 | Tel 855-425-5478 | www.Alliqua.com

NDC#58189-044-01



PRODUCT CODE
AB-400400-001



(01) 00816789010006

HYDROGEL DRESSING

STERILE
SINGLE USE ONLY



Contents: One (1) dressing
4" X 4" (10cm X 10cm)

Hydress™ should be placed on the wound and bandaged as necessary.

Directions: Remove dressing from packaging and remove the blue plastic sheet from dressing. Apply exposed side of dressing to wound and bandage as necessary. Change dressing daily.

For use on minor burns, minor cuts and scrapes.

Ingredients: Deionized water, polyethylene oxide, methylparaben, propylparaben and citric acid.

To facilitate removal, moisten with sterile water.

Use before expiration date.

Manufactured by AquaMed Technologies Inc. for Alliqua BioMedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047 | www.Alliqua.com

Principal Display Panel – 4" X 5" Dressing Box Label

NDC 58198-145-10

SILVER

SEAL®

PRODUCT CODE

AB-400500-105

(01)10816789010010

HYDROGEL DRESSING

STERILE Contents: Ten (10) dressings

SINGLE USE ONLY 4" X 5" (10cm X 12.7cm)

Manufactured by Aquamed Technologies Inc. for Alliqua Biomedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047 | www.Alliqua.com

NDC#58189-145-10



PRODUCT CODE
AB-400500-105



(01)10816789010010

HYDROGEL DRESSING

.....

STERILE
SINGLE USE ONLY

Contents: Ten (10) dressings

4" X 5" (10cm X 12.7cm)

Apply as a single layer.

Overlay with conventional dressings.

SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier.

Applicable to:

partial and full thickness dermal ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.



SilverSeal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber.

Antibacterial activity was confirmed via in vitro studies of bacterial growth in the dressing.

To facilitate removal, moisten with sterile water.

**Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.
Use before expiration date.**

Manufactured by AquaMed Technologies Inc. for Alliqua BioMedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047 | www.Alliqua.com

Principal Display Panel – 4" X 5" Dressing Box Label

NDC 58198-145-01

SILVER
SEAL®

PRODUCT CODE
AB-400500-105

(01)00816789010013

HYDROGEL DRESSING

STERILE Contents: One (1) dressing
SINGLE USE ONLY 4" X 5" (10cm X 12.7cm)


Manufactured by Aquamed Technologies Inc. for Alliqua Biomedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047 | www.Alliqua.com

NDC#58189-145-01



PRODUCT CODE

AB-400500-105



(01) 00816789010013

.....HYDROGEL DRESSING.....

STERILE
SINGLE USE ONLY

Contents: One (1) dressing

4" X 5" (10cm X 12.7cm)

Apply as a single layer.

Overlay with conventional dressings.

SilverSeal® Hydrogel Dressing is intended for the management of wounds and to provide an antimicrobial barrier.

Applicable to:

partial and full thickness dermal ulcers, leg ulcers (vascular, venous, pressure and diabetic), superficial wounds, abrasions, first and second degree burns, and donor sites.

SilverSeal® Hydrogel Dressing contains approximately 1.5% X-Static® silver coated fiber.

Antibacterial activity was confirmed via in vitro studies of bacterial growth in the dressing.

To facilitate removal, moisten with sterile water.



**Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.
Use before expiration date.**

Manufactured by AquaMed Technologies Inc. for Alliqua BioMedical Inc.
2150 Cabot Boulevard West, Langhorne, PA 19047 | www.Alliqua.com

X-STATIC SILVERSEAL HYDROGEL DRESSING 4 X 5

dressing, wound and burn, hydrogel w/drug and/or biologic

Product Information			
Product Type	MEDICAL DEVICE	Item Code (Source)	NHRIC:58 189-145

Packaging				
#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:58 189-145-10	10 in 1 CARTON		
1	NHRIC:58 189-145-01	1 in 1 POUCH		

Marketing Information			
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
Premarket Notification	K040019	05/15/2013	

X-STATIC SILVERSEAL HYDROGEL DRESSING 2 X 3			
dressing, wound and burn, hydrogel w/drug and/or biologic			

Product Information			
Product Type	MEDICAL DEVICE	Item Code (Source)	NHRIC:58 189-123

Packaging				
#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:58 189-123-10	10 in 1 CARTON		
1	NHRIC:58 189-123-01	1 in 1 POUCH		

Marketing Information			
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
Premarket Notification	K040019	05/15/2013	

X-STATIC SILVERSEAL HYDROGEL DRESSING 2 X 2			
dressing, wound and burn, hydrogel w/drug and/or biologic			

Product Information			
Product Type	MEDICAL DEVICE	Item Code (Source)	NHRIC:58 189-122

Packaging				
#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:58 189-122-10	10 in 1 CARTON		
1	NHRIC:58 189-122-01	1 in 1 POUCH		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
Premarket Notification	K040019	06/01/2013	

X-STATIC SILVERSEAL HYDROGEL DRESSING 3

dressings, wound and burn, hydrogel w/drug and/or biologic

Product Information

Product Type	MEDICAL DEVICE	Item Code (Source)	NHRIC:58189-130
--------------	----------------	--------------------	-----------------

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:58189-130-10	10 in 1 CARTON		
1	NHRIC:58189-130-01	1 in 1 POUCH		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
Premarket Notification	K040019	06/01/2013	

HYDRESS

dressings, wound, hydrophilic

Product Information

Product Type	MEDICAL DEVICE	Item Code (Source)	NHRIC:58189-144
--------------	----------------	--------------------	-----------------

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:58189-144-10	10 in 1 CARTON		
1	NHRIC:58189-144-01	1 in 1 POUCH		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
Exempt Device	FRO	05/15/2013	

Labeler - Alliqua Inc (872967679)

