

**QCARE- cetylpyridinium chloride**  
**Sage Products LLC**

*Disclaimer: Most OTC drugs are not reviewed and approved by FDA, however they may be marketed if they comply with applicable regulations and policies. FDA has not evaluated whether this product complies.*

-----

**Corinz™ Antiseptic Cleansing and Moisturizing Oral Rinse**

**Drug Facts**

<b>Active ingredient:</b>	<b>Purpose</b>
CORINZ ANTISEPTIC CLEANSING AND MOISTURIZING ORAL RINSE:	
Cetylpyridinium chloride 0.05%	Antiseptic Rinse

**Uses**

- Aids in the removal of secretions and debris and helps reduce the chance of infection in minor oral irritation.

**Warnings**

**Stop use and ask a doctor if:**

- Sore mouth symptoms do not improve in 7 days.
- Swelling, rash or fever develops.
- Irritation, pain or redness persists or worsens.

**Keep out of reach of children.**

If more than used for antisepsis is accidentally swallowed, get medical help or contact a Poison Control Center right away.

**Directions**

**Toothette Oral Care Suction Toothbrush**

- **Before opening**, turn package over, burst solution packet with thumbs.
- Peel lid to open.
- Attach Toothbrush to suction line.
- Clean teeth and oral cavity for approximately one minute.
- To suction, place thumb over port.
- To clear tubing, rinse with sterile saline or appropriate solution.
- Discard Toothbrush. Reattach Covered Yankauer to suction line.
  
- Use up to 4 times daily or as directed by a dentist or doctor.
- Children under 12 years of age: supervise use.
- Children under 3 years of age: consult a dentist or doctor.
- Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
- Ensure foam is intact after use. If not, remove any particles from oral cavity.

**Toothette Oral Care Suction Swab**

- **Before opening**, turn package over, burst solution packet with thumbs.

- Peel lid to open.
- Attach Suction Swab to suction line.
- Clean teeth and oral cavity for approximately one minute.
- To suction, place thumb over port.
- To clear tubing, rinse with sterile saline or appropriate solution.
- Discard Suction Swab. Reattach Covered Yankauer to suction line.
  
- Use up to 4 times daily or as directed by a dentist or doctor.
- Children under 12 years of age: supervise use.
- Children under 3 years of age: consult a dentist or doctor.
- Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
- Ensure foam is intact after use. If not, remove any particles from oral cavity.

***Inactive ingredients***

Water, glycerin, xylitol, spearmint flavor, potassium sorbate, polysorbate 20, polysorbate 80, hydroxyethylcellulose, citric acid, sodium saccharin, menthol.

***Questions?***

Call toll-free 800-323-2220.

LATEX FREE • FOR SINGLE USE ONLY • MADE IN U.S.A.

Tear Here

# Q-Care®


# Continue care™

for the non-ventilated patient



### Setup Instructions:

- 1** Connect tubing to suction port on canister.
- 2** Connect the Covered Yankauer to tubing.
- 3** Connecting tools: Detach the Covered Yankauer from tubing. Attach cleansing and suctioning tools directly to tubing.



Cleansing and Suctioning Tools attach directly to standard suction lines.

Using the Covered Yankauer to aid with removal of secretions:

- Follow hospital protocol.
- Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
- Retract sleeve until cap locks in place.
- To suction, move switch to ON. When finished, move switch to OFF to prevent loss of suction power.
- Rinse after use until tubing is clear.
- Between uses, pull sleeve over tip to cover.
- Discard Yankauer after using for 24 hours, or according to hospital protocol.
- For single patient oral use.

#### Contents:

- 1 Covered Yankauer
- 3 Suction Toothbrush packages with Corinz™ Antiseptic Cleansing and Moisturizing Oral Rinse

#### Ingredients:

Corinz™ Antiseptic Cleansing and Moisturizing Oral Rinse: Cetylpyridinium chloride 0.05%  
Kit also contains sodium bicarbonate. For detailed listing of ingredients, see the enclosed individual packages.

#### Caution:

Federal (U.S.A.) law restricts this device (kit) to sale by or on the order of a physician or licensed practitioner.

Non-sterile • Latex free

## Reorder # 6308



MADE IN U.S.A.  
Patents: [www.sageproducts.com/patents](http://www.sageproducts.com/patents)

# SAGE®

PRODUCTS

Sage Products LLC  
3909 Three Oaks Road | Cary, Illinois 60013  
[www.sageproducts.com](http://www.sageproducts.com) | 800-323-2220

# Q-Care®

# Continue care™

for the non-ventilated patient

featuring  
**Corinz™**  
Antiseptic Cleansing &  
Moisturizing Oral Rinse

## Setup Instructions:


**1** Connect tubing to suction port on canister.

**2** Connect the Covered Yankauer to tubing.

**3** Connecting tools: Detach the Covered Yankauer from tubing. Attach cleansing and suctioning tools directly to tubing.

Cleansing and Suctioning Tools attach directly to standard suction lines.

Using the Covered Yankauer to aid with removal of secretions:



- Follow hospital protocol.
- Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
- Retract sleeve until cap locks in place.
- To suction, move switch to ON. When finished, move switch to OFF to prevent loss of suction power.
- Rinse after use until tubing is clear.
- Between uses, pull sleeve over tip to cover.
- Discard Yankauer after using for 24 hours, or according to hospital protocol.
- For single patient oral use.

### Contents:

- 1 Covered Yankauer
- 2 Suction Toothbrush packages with Corinz™ Antiseptic Cleansing and Moisturizing Oral Rinse
- 2 Suction Oral Swab packages with Corinz™ Antiseptic Cleansing and Moisturizing Oral Rinse

### Ingredients:

Corinz™ Antiseptic Cleansing and Moisturizing Oral Rinse: Cetylpyridinium chloride 0.05%  
Kit also contains sodium bicarbonate. For detailed listing of ingredients, see the enclosed individual packages.

### Caution:

Federal (U.S.A.) law restricts this device (kit) to sale by or on the order of a physician or licensed practitioner.

Non-sterile • Latex free

Reorder # 6306



MADE IN U.S.A.  
Patents: [www.sageproducts.com/patents](http://www.sageproducts.com/patents)

# SAGE®

PRODUCTS

Sage Products LLC  
3909 Three Oaks Road | Cary, Illinois 60013  
[www.sageproducts.com](http://www.sageproducts.com) | 800-323-2220



# Suction Toothbrush

with Corinz™ Antiseptic Cleansing and Moisturizing Oral Rinse

CONTENTS: 1 Suction Toothbrush with sodium bicarbonate, 0.25 fl. oz./7ml Corinz Antiseptic Cleansing and Moisturizing Oral Rinse

## Drug Facts

Active ingredient	Purpose
CORINZ ANTISEPTIC CLEANSING AND MOISTURIZING ORAL RINSE: Cetylpyridinium chloride 0.05%.....	Antiseptic Rinse

## Uses

- Aids in the removal of secretions and debris and helps reduce the chance of infection in minor oral irritation.

## Warnings

Stop use and ask a doctor if:

- Sore mouth symptoms do not improve in 7 days.
- Swelling, rash or fever develops.
- Irritation, pain or redness persists or worsens.

Keep out of reach of children.

If more than used for antisepsis is accidentally swallowed, get medical help or contact a Poison Control Center right away.

## Directions

- **Before opening**, turn package over, burst solution packet with thumbs.
- Peel lid to open.
- Attach Toothbrush to suction line.
- Clean teeth and oral cavity for approximately one minute.
- To suction, place thumb over port.
- To clear tubing, rinse with sterile saline or appropriate solution.
- Discard Toothbrush. Reattach Covered Yankauer to suction line.

- Use up to 4 times daily or as directed by a dentist or doctor.
- Children under 12 years of age: supervise use.
- Children under 3 years of age: consult a dentist or doctor.
- Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
- Ensure foam is intact after use. If not, remove any particles from oral cavity.

## Inactive ingredients

Water, glycerin, xylitol, spearmint flavor, potassium sorbate, polysorbate 20, polysorbate 80, hydroxyethylcellulose, citric acid, sodium saccharin, menthol.

Questions? Call toll-free 800-323-2220.



# Suction Swab

with Corinz™ Antiseptic Cleansing and Moisturizing Oral Rinse

CONTENTS: 1 Suction Swab with sodium bicarbonate, 0.25 fl. oz./7ml Corinz Antiseptic Cleansing and Moisturizing Oral Rinse

## Drug Facts

Active ingredient	Purpose
CORINZ ANTISEPTIC CLEANSING AND MOISTURIZING ORAL RINSE: Cetylpyridinium chloride 0.05%.....	Antiseptic Rinse

## Uses

- Aids in the removal of secretions and debris and helps reduce the chance of infection in minor oral irritation.

## Warnings

Stop use and ask a doctor if:

- Sore mouth symptoms do not improve in 7 days.
- Swelling, rash or fever develops.
- Irritation, pain or redness persists or worsens.

Keep out of reach of children.

If more than used for antisepsis is accidentally swallowed, get medical help or contact a Poison Control Center right away.

## Directions

- **Before opening**, turn package over, burst solution packet with thumbs.
- Peel lid to open.
- Attach Suction Swab to suction line.
- Clean teeth and oral cavity for approximately one minute.
- To suction, place thumb over port.
- To clear tubing, rinse with sterile saline or appropriate solution.
- Discard Suction Swab. Reattach Covered Yankauer to suction line.

- Use up to 4 times daily or as directed by a dentist or doctor.
- Children under 12 years of age: supervise use.
- Children under 3 years of age: consult a dentist or doctor.
- Use a bite block when performing oral care on patients with altered levels of consciousness or those who cannot comprehend commands.
- Ensure foam is intact after use. If not, remove any particles from oral cavity.

## Inactive ingredients

Water, glycerin, xylitol, spearmint flavor, potassium sorbate, polysorbate 20, polysorbate 80, hydroxyethylcellulose, citric acid, sodium saccharin, menthol.

Questions? Call toll-free 800-323-2220.

**Sodium Bicarbonate ingredients:** Water, sodium bicarbonate, sodium carboxymethyl cellulose, sodium lauryl sulfate, flavor, sodium saccharin, sodium benzoate.

LATEX FREE • FOR SINGLE USE ONLY • MADE IN U.S.A.



**Oral Check™** | Scan bar code for compliance to protocol

**Sodium Bicarbonate ingredients:** Water, sodium bicarbonate, sodium carboxymethyl cellulose, sodium lauryl sulfate, flavor, sodium saccharin, sodium benzoate.

LATEX FREE • FOR SINGLE USE ONLY • MADE IN U.S.A.



**Oral Check™** | Scan bar code for compliance to protocol



**Antiseptic Cleansing &  
Moisturizing Oral Rinse**

**Active Ingredient:  
Cetylpyridinium  
chloride 0.05%**



**BURST PACKET  
BEFORE OPENING  
OUTER PACKAGE**

**Do not remove  
from package  
0.25 fl. oz./7 mL**

R5439



**Antiseptic Cleansing &  
Moisturizing Oral Rinse**

**Active Ingredient:  
Cetylpyridinium  
chloride 0.05%**



**BURST PACKET  
BEFORE OPENING  
OUTER PACKAGE**

**Do not remove  
from package  
0.25 fl. oz./7 mL**



**107**

**QCARE**

cetylpyridinium chloride kit

**Product Information**

<b>Product Type</b>	HUMAN OTC DRUG	<b>Item Code (Source)</b>	NDC:53462-306
---------------------	----------------	---------------------------	---------------

**Packaging**

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:53462-306-60	1 in 1 KIT; Type 4: Device Coated/Impregnated/Otherwise Combined with Drug		

**Quantity of Parts**

Part #	Package Quantity	Total Product Quantity
Part 1	4 POUCH	28 mL

**Part 1 of 1****CORINZ**

cetylpyridinium chloride rinse

**Product Information**

<b>Item Code (Source)</b>	NDC:53462-375
<b>Route of Administration</b>	ORAL

**Active Ingredient/Active Moiety**

Ingredient Name	Basis of Strength	Strength
CETYLPYRIDINIUM CHLORIDE (UNII: D9OM4SK49P) (CETYLPYRIDINIUM - UNII:CUB7J10JV3)	CETYLPYRIDINIUM CHLORIDE	.5 mg in 1 mL

**Inactive Ingredients**

Ingredient Name	Strength
WATER (UNII: 059QF0KO0R)	
GLYCERIN (UNII: PDC6A3C0OX)	
XYLITOL (UNII: VCQ006KQ1E)	
POTASSIUM SORBATE (UNII: 1VPU26JZZ4)	
POLYSORBATE 20 (UNII: 7T1F30V5YH)	
POLYSORBATE 80 (UNII: 6OZP39ZG8H)	
HYDROXYETHYL CELLULOSE (4000 MPAS AT 1%) (UNII: ZYD53NBL45)	
CITRIC ACID MONOHYDRATE (UNII: 2968PHW8QP)	
SACCHARIN SODIUM (UNII: SB8ZUX40TY)	
MENTHOL, UNSPECIFIED FORM (UNII: L7T10EIP3A)	



**Packaging**

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1		1 in 1 PACKET		
1		7 mL in 1 POUCH		

**Marketing Information**

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
OTC monograph not final	part356	12/28/2001	

**Marketing Information**

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
OTC monograph not final	part356	12/28/2001	

**QCARE**

cetylpyridinium chloride kit

**Product Information**

Product Type	Item Code (Source)
HUMAN OTC DRUG	NDC:53462-308

**Packaging**

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:53462-308-60	1 in 1 KIT; Type 4: Device Coated/Impregnated/Otherwise Combined with Drug		

**Quantity of Parts**

Part #	Package Quantity	Total Product Quantity
Part 1	3 POUCH	21 mL

**Part 1 of 1****CORINZ**

cetylpyridinium chloride rinse

**Product Information**

Item Code (Source)	NDC:53462-375
Route of Administration	ORAL

Active Ingredient/Active Moiety		
Ingredient Name	Basis of Strength	Strength
CETYLPYRIDINIUM CHLORIDE (UNII: D9OM4SK49P) (CETYLPYRIDINIUM - UNII:CUB7JJ0JV3)	CETYLPYRIDINIUM CHLORIDE	.5 mg in 1 mL

Inactive Ingredients	
Ingredient Name	Strength
WATER (UNII: 059QF0KO0R)	
GLYCERIN (UNII: PDC6A3C0OX)	
XYLITOL (UNII: VCQ006KQ1E)	
POTASSIUM SORBATE (UNII: 1VPU26JZZ4)	
POLYSORBATE 20 (UNII: 7T1F30V5YH)	
POLYSORBATE 80 (UNII: 6OZP39ZG8H)	
HYDROXYETHYL CELLULOSE (4000 MPAS AT 1%) (UNII: ZYD53NBL45)	
CITRIC ACID MONOHYDRATE (UNII: 2968PHW8QP)	
SACCHARIN SODIUM (UNII: SB8ZUX40TY)	
MENTHOL, UNSPECIFIED FORM (UNII: L7T10EIP3A)	

Packaging				
#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1		1 in 1 PACKET		
1		7 mL in 1 POUCH		

Marketing Information			
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
OTC monograph not final	part356	12/28/2001	

Marketing Information			
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
OTC monograph not final	part356	12/28/2001	

**Labeler** - Sage Products LLC (054326178)

**Registrant** - Sage Products LLC (054326178)

Establishment			
Name	Address	ID/FEI	Business Operations
Sage Products LLC		054326178	MANUFACTURE(53462-308, 53462-306)