


NDA 208799/S-026
 Page 3

NDA 208799/S-26 dated May 24, 2023: 55 mcg/14 mcg




NDA 208799/S-026
Page 4




NDA 208799/S-026
Page 6

Get to know your AirDuo® Dighaler® dry powder inhaler.
The **yellow cap** on your inhaler is used to load a dose. Every time the yellow cap is opened and it clicks, one dose is ready to be inhaled.
Read the full Instructions for Use before you start.


Open



Inhale



Close



Important:
Always close the yellow cap to re-load your next dose.

Open the yellow cap fully until you feel and hear a click.


Breathe in quickly and deeply to release medication and hold your breath for about 10 seconds.

Close the yellow cap to re-load your next dose.


See more info below

A built-in electronic module records and stores information about your inhaler events that you can view in the app. AirDuo Dighaler does not need to be connected to the app in order for you to take your medicine.

Scan QR code for each new Dighaler to see your Inhaler events.



Use the app to view and share your Dighaler event data.



For information on your Dighaler and app, call 1-877-410-DIGI (3444) or go to www.dighaler.com.

NDA 208799/S-026
 Page 7

NDA 208799/S-26 dated May 24, 2023: 232 mcg/14 mcg

Non-Varnish Area For Serialization Coating Area

NDC 59310-136-06

AirDuo[®] digihaler[®] 232/14
 (fluticasone propionate 232 mcg and salmeterol 14 mcg) Inhalation Powder



AirDuo[®] digihaler[®] 232/14
 (fluticasone propionate 232 mcg and salmeterol 14 mcg)
 Inhalation Powder

CONTENTS: Each metered dose contains 232 mcg of fluticasone propionate and 14 mcg of salmeterol base, equivalent to 20.9 mcg of salmeterol xinafoate, with lactose.

Usual Dosage: See package insert for full prescribing information.

IMPORTANT SAFETY INFORMATION

Do not increase the prescribed dose or discontinue abruptly. While taking AirDuo[®] Digihaler[®], other inhaled medicines should be used only as prescribed by your physician.

NEVER WASH OR PUT ANY PART OF YOUR INHALER IN WATER.

Store inhaler at controlled room temperature between 15°C and 25°C (59°F and 77°F). Avoid exposure to extreme heat, cold, or humidity.

IMPORTANT

Discard inhaler when the dose counter displays 0, 30 days after removal from the foil pouch or after the expiration date on the product, whichever comes first.

For information on your Digihaler and app, see the inside of the carton, call 1-877-410-DIGI (3444), or go to www.digihaler.com.

Manufactured for:
 Teva Pharmaceuticals USA, Inc., Parsippany, NJ 07054

MANUFACTURED IN IRELAND
 ©2023 Teva Respiratory, LLC
 FCC ID: 2A/VSPDR000401

For California Only: This product uses batteries containing Perchlorate. Material - special handling may apply. See www.dtsc.ca.gov/hazardouswaste/perchlorate [22 CCR 67384.4(a)].

OTR 033823 TO US0811



N 59310-136-06 1

NDC 59310-136-06
 Rx Only

AirDuo[®] digihaler[®] 232/14
 (fluticasone propionate 232 mcg and salmeterol 14 mcg)
 Inhalation Powder

Scan the QR code on top of your inhaler to download the Digihaler app.

AirDuo[®] digihaler[®] 232/14
 (fluticasone propionate 232 mcg and salmeterol 14 mcg)
 Inhalation Powder

232 mcg/14 mcg

60 metered inhalations
 0.45 g net contents

Refer to enclosed Patient Leaflet and Instructions for Use for detailed information on product use and handling.

See inside for more information on your AirDuo Digihaler and downloading the app.

OPEN HERE

NDC 59310-136-06
 Rx Only

AirDuo[®] digihaler[®] 232/14
 (fluticasone propionate 232 mcg and salmeterol 14 mcg)
 Inhalation Powder

Scan the QR code on top of your inhaler to download the Digihaler app.

AirDuo[®] digihaler[®] 232/14
 (fluticasone propionate 232 mcg and salmeterol 14 mcg)
 Inhalation Powder

232 mcg/14 mcg

60 metered inhalations
 0.45 g net contents

Refer to enclosed Patient Leaflet and Instructions for Use for detailed information on product use and handling.

See inside for more information on your AirDuo Digihaler and downloading the app.

OPEN HERE


NDA 208799/S-026
Page 8

Get to know your AirDuo[®] DigiInhaler[™] dry powder inhaler.
The **yellow cap** on your inhaler is used to load a dose. Every time the yellow cap is opened and it clicks, one dose is ready to be inhaled.

Read the full Instructions for Use before you start.


Important:
Always close the yellow cap to re-load your next dose.

Open




CLICK

Inhale



Close




Breathe in quickly and deeply to release medication and hold your breath for about 10 seconds.

Open the yellow cap fully until you feel and hear a click.


Close the yellow cap to re-load your next dose.

See more info below

Scan QR code for each new DigiInhaler to see your inhaler events.



Use the app to view and share your DigiInhaler event data.



A built-in electronic module records and stores information about your inhaler events that you can view in the app. AirDuo DigiInhaler does not need to be connected to the app in order for you to take your medicine.

For information on your DigiInhaler and app, call 1-877-410-DIGI (3444) or go to www.digiinhaler.com.

This label may not be the latest approved by FDA.
For current labeling information, please visit <https://www.fda.gov/drugsatfda>



David
Lewis

Digitally signed by David Lewis
Date: 11/16/2023 05:12:45PM
GUID: 508da72000029f287fa31e664741b577