

DEZINFEKCE DISINFECTANT- hypochlorous acid disinfectant liquid
Shanghai JiYi Biological Technology Co., Ltd.

Disclaimer: This drug has not been found by FDA to be safe and effective, and this labeling has not been approved by FDA. For further information about unapproved drugs, [click here](#).

JiYi, DEZINFEKCE Disinfectant, Liquid, 25L

Active Ingredient(s)

Hypochlorous acid 700 mg/L. Purpose: Kill microorganism.

Purpose

The product can kill common bacteria in hospital infection, pathogenic entero bacteria, pyogenic coccus, pathogenic yeast, escherichia coli, staphylococcus aureus, pseudomonas aeruginosa, candida albicans, bacillus subtilis var. niger spores and other microorganisms.

Use

Use the product to spray, wipe, soak for disinfection. Air dry is recommended and here is no need for secondary cleaning.

Scope of application:

Surface disinfection: spray and wipe to disinfect stair handrails, door handles, desktop, seats, floors, etc.;

Disinfection of fruits and vegetables: soak fruits and vegetables to disinfect and kill surface microorganisms.

Warnings

External Use Only

- Do not replace containers
- Do not mix it with other disinfection products
- Keep it sealed, and use it as soon as possible after unpacking
- Spraying on the fabrics is not recommended
- Not suitable for articles made of copper or carbon steel
- Keep out of reach children

Do not use

- on children less than 2 months of age
- on open skin wounds
- mixing this product with other disinfectants

When using this product keep out of eyes, ears, and mouth.

In case of contact with eyes, rinse eyes thoroughly with water.

This product is for external use. It should not be taken orally.

This product has a certain corrosive effect on metal objects other than stainless steel, and should be used with caution.

When sterilizing articles containing more organic objects, they should be completely removed and then disinfected.

Use with caution for those with sensitive skin.

Do not mix this product with other disinfectants.

Stop use and ask a doctor if irritation or rash occurs. These may be signs of a serious condition.

Directions

1. Surface disinfection of hard objects: Use the stock solution to wipe and disinfect, act for 20 minutes, and air-dry naturally without secondary cleaning.
2. Disinfection of food and drinking utensils: Soak it in the product for 30 minutes and then take it out, and air-dry naturally without secondary cleaning.
3. Air disinfection: Use the stock solution to atomize the indoor air at an amount of 10mL/m³ for 30 minutes.


Other information

Storage: place it in a normal temperature and light-shielded environment.

Inactive ingredients

Sodium chloride, Water

Package Label - Principal Display Panel



DEZINFEKCE™ DISINFECTANT



Authorized by AE Company, German
100% German Technologies

Please read the instructions carefully
and use it according to the instructions

Product name: DEZINFEKCE™ Disinfectant

Active Ingredient and content: Hypochlorous acid

Category of microorganism killed: This product can kill common bacteria in hospital infection, pathogenic entero bacteria, pyogenic coccus, pathogenic yeast, escherichia coli, staphylococcus aureus, pseudomonas aeruginosa, candida albicans, bacillus subtilis var. niger spores and other microorganisms.

Usage: Use the product to spray, wipe, soak for disinfection. Air dry is recommended and here is no need for secondary cleaning.

Scope of application:

- Surface disinfection: spray and wipe to disinfect stair handrails, door handles, desktop, seats, floors, etc.;
- Disinfection of fruits and vegetables: soak fruits and vegetables to disinfect and kill surface microorganisms.

Notes:

- Do not replace containers
- Do not mix it with other disinfection products
- Keep it sealed, and use it as soon as possible after unpacking
- Spraying on the fabrics is not recommended
- Not suitable for articles made of copper or carbon steel
- Keep out of reach children

Product features: The product is non-toxic through the mouth and has no irritation to skin.

Dosage form: Liquid

Specifications: 25L/barrel

Storage: Place it in a normal temperature and light-shielded environment.

Validity: 12 months

Production date: See the barrel

Production batch No.: See the barrel

Executive standard: Q/SHJY003

Approval No.: SWX ZI (2018) No. 3208-0001

Responsible: Shanghai JIYI Biotech Co., Ltd.

Address: Room 5700, Building 1, No. 1150, Lanfeng Road, Fengxian District, Shanghai City, China

Postal code: 200082

Tel.: 086-021-80198556

Service hotline: 086-400-928-2766

Qualified
Inspected



*Direct spraying



* Pour it into humidifier



JIYI BIOTECH



AQUA
ELECTRA

Made in China

DEZINFEKCE DISINFECTANT

hypochlorous acid disinfectant liquid

Product Information

Product Type	HUMAN OTC DRUG	Item Code (Source)	NDC:75316-125
Route of Administration	EXTRACORPOREAL		

Active Ingredient/Active Moiety				
Ingredient Name			Basis of Strength	Strength
HYPOCHLOROUS ACID (UNII: 712K4CDC10) (HYPOCHLOROUS ACID - UNII:712K4CDC10)			HYPOCHLOROUS ACID	17.5 g in 25 L
Inactive Ingredients				
Ingredient Name			Strength	
WATER (UNII: 059QF0KO0R)				
SODIUM CHLORIDE (UNII: 451W47IQ8X)				
Packaging				
#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:75316-125-16	32 in 1 TRAY	04/23/2020	
1	NDC:75316-125-01	25 L in 1 DRUM; Type 0: Not a Combination Product		
Marketing Information				
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date	
unapproved drug other		04/23/2020		

Labeler - Shanghai Jiyi Biological Technology Co., Ltd. (554529167)

Registrant - A03 Lab of BTS (548009541)

Establishment

Name	Address	ID/FEI	Business Operations
Shanghai Jiyi Biological Technology Co., Ltd.		554529167	label(75316-125) , manufacture(75316-125)

Revised: 4/2020

Shanghai Jiyi Biological Technology Co., Ltd.