

VIAL KIT3.98- chitosan liquid Oxytical

Disclaimer: This drug has not been found by FDA to be safe and effective, and this labeling has not been approved by FDA. For further information about unapproved drugs, click here.

Active Ingredients

Chitosan

Purpose

Skin regeneration and Wrinkle relief

Use

The skin's surface appears firmer and smoother with a more youthful look. Remodels the facial contours and reduces the thickness of subcutaneous fatty tissue helping to reduce the "double chin effect" providing you with a youthful complexion. Reinforces skin elasticity to provide plump and radiant complexion.

Directions

[Preparation stage]

After cleansing, use toner to clean up the skin.

It must be used only after cleansing is complete, and especially pay attention to the cleanliness of your hands.

If you use it after exfoliation, you can get more satisfactory results, and it is more effective if you give a simple massage before use.

[Application stage]

1. Prepare Core peptide(1st), Stable Chitin(2nd) and a syringe.
2. Core Peptide(1st) only opens the plastic cap.
3. A syringe inlet is injected into the center of the hole at the top of the Core Peptide(1st) to suction. (Total 3ml based on syringe).

After the syringe is removed, the inlet is automatically sealed to prevent outside air from entering.

4. Open the Stable Chitin(2nd) after referring to the precautions and inject the Core Peptide(1st).

Due to the nature of the Stable Chitin component, it is not soluble in ampoule other than Core Peptide. (Not compatible)

5. Close the vial opening of the Stable Chitin(2nd) with an elastomer tube and shake to dissolve the Stable Chitin.

Check with eyes and if the Stable Chitin component does not appear in the form of a thread (slightly thick), apply it to the skin.

When Using

Although it can be used immediately after melting, it is recommended to use after 12 hours after mixing the product and closing it with an elastomer tube after mixing the product to increase the activity and maximize the efficacy due to the nature of the formulation.

After blending, only the elastomer tube is closed and stored at room temperature fresh out of direct sunlight.

(Up to 21 days)

Ingredient

Chitosan

Warnings

For external use only.

Keep away from heat.

Stop Using

Stop use and ask a physician if skin irritation and redness occurs more than 72 hours.

Keep out of reach of children

Keep out of reach of children, if swallowed, get medical help right away.

Package Label

Oxytical Vial Kit 3.98

Normalizes skin balance / Your skin becomes supple, soft and resilient with a radiant finish.

Patent application No. : 10-2019-0104521

Ingredients active Saussurea involucreta help to improve the appearance of firmness and elasticity on the skin's surface. Bee Venom combined with high-grade Peptides helps to visibly reduce the appearance of fine lines and wrinkles on the skin's surface. Sodium Hyaluronate, known for its restorative and anti-aging abilities, revitalizes and moisturizes skin's outer surface layers to improve appearance of fine lines and wrinkles.

Use The skin's surface appears firmer and smoother with a more youthful look. Remodels the facial contours and reduces the thickness of subcutaneous fatty tissue helping to reduce the "double chin effect" providing you with a youthful complexion. Reinforces skin elasticity to provide plump and radiant complexion.

Warnings

1) For external use only. Keep away from heat 2) Stop use and ask a physician if skin irritation and redness occurs more than 72 hours. 3) Keep out of reach of children, if swallowed, get medical help right away.

How to use : Picture reference

Direction

Although it can be used immediately after melting, it is recommended to use after 12 hours after mixing the product and closing it with an elastomer tube after mixing the product to increase the activity and maximize the efficacy due to the nature of the formulation. After blending, only the elastomer tube is closed and stored at room temperature fresh out of direct sunlight.(Up to 21 days)

Ingredients

1. chitosan 2. Water, Butylene Glycol, Propanediol, Glycereth-26, Saussurea Involucrata Extract, Sodium Hyaluronate, Acetyl Hexapeptide-8, Copper Tripeptide-1, Sodium Lactate, PEG-60 Hydrogenated Castor Oil, Chlorophenesin, Ethylhexylglycerin, Lactic Acid, Fragrance, Dipeptide Diaminobutyryl Benzylamide Diacetate, Bee Venom, CI 19140

Customer care center : 070-5225-5990 **Manufacturer** : ENBIOSCIENCE.,LTD 33, Dodam 4-ro, Seo-gu, Incheon, Republic of Korea **Distributor** : Oxytical 233ho, Maeheon-ro-46, Seocho-gu, Seoul, Republic of Korea

Manufacturing number and date of expiration : Notified separately

* This product can be exchanged or compensated according to notification by the Fair Trade Commission.

VIAL KIT3.98

chitosan liquid

Product Information

Product Type	HUMAN OTC DRUG	Item Code (Source)	NDC:80851-302
Route of Administration	TOPICAL		

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
CHITOSAN MEDIUM MOLECULAR WEIGHT (200-400 MPA.S) (UNII: 5GV09YMO52) (CHITOSAN MEDIUM MOLECULAR WEIGHT (200-400 MPA.S) - UNII:5GV09YMO52)	CHITOSAN MEDIUM MOLECULAR WEIGHT (200-400 MPA.S)	99.99 mg in 100 mL

Inactive Ingredients

Ingredient Name	Strength
WATER (UNII: 059QF0KO0R)	0.001 mL in 100 mL

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:80851-302-01	0.02 mL in 1 BOTTLE, GLASS; Type 0: Not a Combination Product	11/12/2020	

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
unapproved drug other		11/12/2020	

Labeler - Oxytical (695504108)

Registrant - Oxytical (695504108)

Establishment

Name	Address	ID/FEI	Business Operations
Oxytical		695504108	manufacture(80851-302)

Revised: 11/2022

Oxytical