TRICHLOROISOCYANURIC ACID- trichlorois o cyanuric acid liquid GUANGZHOU YIQI TECHNOLOGY CO.,LTD

Disclaimer: This drug has not been found by FDA to be safe and effective, and this labeling has not been approved by FDA. For further information about unapproved drugs, click here.

78295-010 78295-009 replaced by 78295-010 Trichlorois ocyanuric Acid

Active ingredient

Active ingredient: Chlorine 91%.

Purpose

Purpose: Antimicrobial.

WARNINGS

[WARNINGS]

Corrosive I. Store in the orginal Container. For external use only.

Usage:Suitable for air disinfecting.household, public place.Surface of substance, such as white fabrics, tableware

Application	Chlorine(mg/L)	Dilution Rate	method	Time
Tableware	900	1:20	Soak, water flushing	20 Minute
Home Public	300	1:60	Wipe then rinse	20 Minute
Fabric	300	1:60	wipc then rinse	20 Minute

/

When using this product

[When using this product]

Do not use on sikl wool nylon leather paint surfaces Do not use for carbon steel products aluminum copper suitable for disinfection of stainless steel appliances. This product can bleach fabric to a certain extent and only be used for white fabric Do not mix with acid cleaning products and other chemicals. If you splash your eyes or skin rinse with water immediately. If necessary please see a doctor. After disinfecting the tableware please rinse with clean water. Keep out of reach for children.

/

Directions

[Directions]

We suggest to wear gloves when using it To avoid touch directly. keep away from sunlight and heat lain.

Other information

[Other information] Store below 106°F.(41C) May discolor certain fabrics or surfaces.

Inactive ingredient

Inactive ingredient: water

Package Label - Principal Display Panel

78295-010-01 25KGS 1 BAG



Product name: Trichloroisocyanuri Acid **Active ingredient:** Chlorine 91%.

Purpose: Antimicrobial.

Usage: Suitable for air disinfecting, household, public place. Surface of substance, such as white fabrics,tableware

USAGE

Application	Chlorine (mg/L)	Dilution Rate	Method	Time
Tableware	900	1:20	Soak, water flushing	20 Minute
Home Public	300	1:60	Wipe then rinse	20 Minute
Fabric	300	1:60	Wipe then Rinse	20 Minute

TIPS

[WARNINGS] Corrosive, Store in the original

Container . For external use only. [When using this product]Do not use on silk, wool,nylon,leather,paint surfaces. Do not use for carbon steel products, aluminum, copper suitable for and other chemicals. If you splash your eyes or skin, rinse with water immediately. If necessary please see with clean water. Keep out of reach for children. [Directions] We suggest to wear gloves when

disinfection of stainless steel appliances, This product can bleach fabric to a certain extent and only be used for white fabric, Do not mix with acidcleaning products a doctor. After disinfecting the tableware, please rinse

using it, To avoid touch directly. keep away from sunlight and heat, rain,

Trichloroisocyanuric Acid

[Other information] Store below 106°F.(41°C) May discolor certain fabrics or surfaces.

Inactive ingredient: water

[Manufacturer] Guangzhou YIQI Technology

[Address] 10 Lixin Twelve Road, Xintang Town, Zengcheng District, Guangzhou City, China **Quality Guarantee 12 Months.**

MFG:02/04/2020 EXP:01/04/2021











MADE IN CHINA

78295-010-02 50KGS 1 DRUM



Drug Facts

Product name: Trichloroisocyanuri Acid Active ingredient: Chlorine 91%.

Purpose: Antimicrobial.

Usage: Suitable for air disinfecting, household, public place. Surface of substance, such as white fabrics,tableware

USAGE

Application	Chlorine (mg/L)	Dilution Rate	Method	Time
Tableware	900	1:20	Soak, water flushing	20 Minute
Home Public	300	1:60	Wipe then rinse	20 Minute
Fabric	300	1:60	Wipe then Rinse	20 Minute

TIPS

[WARNINGS] Corrosive, Store in the original Container . For external use only.

[When using this product]Do not use on silk, wool,nylon,leather,paint surfaces, Do not use for carbon steel products, aluminum, copper suitable for disinfection of stainless steel appliances, This product can bleach fabric to a certain extent and only be used for white fabric, Do not mix with acidcleaning products and other chemicals. If you splash your eyes or skin, rinse with water immediately. If necessary please see

Trichloroisocyanuric Acid

a doctor. After disinfecting the tableware, please rinse with clean water. Keep out of reach for children. [Directions] We suggest to wear gloves when using it, To avoid touch directly. keep away from sunlight and heat, rain,

[Other information] Store below 106°F.(41°C) May discolor certain fabrics or surfaces.

Inactive ingredient: water

[Manufacturer] Guangzhou YIQI Technology Co., Ltd
[Address] 10 Lixin Twelve Road, Xintang Town, Zengcheng District, Guangzhou City, China Quality Guarantee 12 Months.

MFG:02/04/2020 EXP:01/04/2021

99.99%
STERILIZINGRATE

STERILIZINGRATE

SOKGS/DRUM

MADE IN CHINA

78295-010-03 100KGS 1 DRUM



wool,nylon,leather,paint surfaces. Do not use for

Trichloroisocyanuric Acid

disinfection of stainless steel appliances, This product can bleach fabric to a certain extent and only be used for white fabric, Do not mix with acidcleaning products and other chemicals. If you splash your eyes or skin, rinse with water immediately. If necessary please see a doctor. After disinfecting the tableware, please rinse with clean water. Keep out of reach for children. [Directions] We suggest to wear gloves when using it, To avoid touch directly. keep away from sunlight and heat, rain,

[Other information] Store below 106°F.(41°C) May discolor certain fabrics or surfaces.

Inactive ingredient: water

[Manufacturer] Guangzhou YIQI Technology

[Address] 10 Lixin Twelve Road, Xintang Town, Zengcheng District, Guangzhou City, China Quality Guarantee 12 Months.

MFG:02/04/2020 EXP:01/04/2021



99.99% **STERILIZINGRATE**







MADE IN CHINA

78295-010-04 1000KGS 1 BAG



Drug Facts

Product name: Trichloroisocyanuri Acid

Active ingredient: Chlorine 91%.

Purpose: Antimicrobial.

Usage: Suitable for air disinfecting, household, public place. Surface of substance, such as white

fabrics,tableware

USAGE

Application	(mg/L)	Dilution Rate	Method	Time
Tableware	900	1:20	Soak, water flushing	20 Minute
Home Public	300	1:60	Wipe then rinse	20 Minute
Fabric	300	1:60	Wipe then Rinse	20 Minute

TIPS

Trichloroisocyanuric Acid

[WARNINGS] Corrosive, Store in the original Container . For external use only.

[When using this product]Do not use on silk, wool,nylon,leather,paint surfaces, Do not use for carbon steel products,aluminum,copper suitable for disinfection of stainless steel appliances,This product can bleach fabric to a certain extent and only be used for white fabric,Do not mix with acidcleaning products and other chemicals. If you splash your eyes or skin, rinse with water immediately. If necessary please see a doctor. After disinfecting the tableware, please rinse with clean water. Keep out of reach for children.

[Directions] We suggest to wear gloves when using it, To avoid touch directly. keep away from sunlight and heat, rain,

[Other information] Store below 106°F.(41°C) May discolor certain fabrics or surfaces.

Inactive ingredient: water

[Manufacturer] Guangzhou YIQI Technology Co., Ltd

[Address] 10 Lixin Twelve Road, Xintang Town, Zengcheng District, Guangzhou City, China Quality Guarantee 12 Months. MFG:02/04/2020 EXP:01/04/2021



99.99% STERILIZINGRATE





TRICHLOROISOCYANURIC ACID

trichloroisocyanuric acid liquid

Product Information

Product Type HUMAN OTC DRUG Item Code (Source) NDC:78295-010

Route of Administration TOPICAL

Active Ingredient/Active Moiety

Ingredient NameBasis of StrengthStrengthCHLORINE (UNII: 4R7X102820) (CHLORINE - UNII:4R7X102820)CHLORINE91 kg in 100 kg

Inactive Ingredients

Ingredient Name Strength

WATER (UNII: 059QF0KO0R)

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:78295-010-01	25 kg in 1 BAG; Type 0: Not a Combination Product	12/15/2020	
2	NDC:78295-010-04	1000 kg in 1 BAG; Type 0: Not a Combination Product	12/15/2020	
3	NDC:78295-010-02	50 kg in 1 DRUM; Type 0: Not a Combination Product	12/15/2020	
4	NDC:78295-010-03	100 kg in 1 DRUM; Type 0: Not a Combination Product	12/15/2020	

Marketing Information				
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date	
unapproved drug other		12/15/2020		

Labeler - GUANGZHOU YIQI TECHNOLOGY CO.,LTD (544707398)

Registrant - GUANGZHOU YIQITECHNOLOGY CO.,LTD (544707398)

Establishment				
Name	Address	ID/FEI	Business Operations	
GUANGZHOU YIQI TECHNOLOGY CO.,LTD		544707398	manufacture(78295-010)	

Revised: 12/2020 GUANGZHOU YIQI TECHNOLOGY CO.,LTD