

TROPHAMINE- isoleucine, leucine, lysine acetate, methionine, phenylalanine, threonine, tryptophan, valine, cysteine hydrochloride, histidine, tyrosine, n-acetyl-tyrosine, alanine, arginine, proline, serine, glycine, aspartic acid, glutamic acid, and taurine solution
B. Braun Medical Inc.

TrophAmine®
(Amino Acid Injection)

**PHARMACY BULK PACKAGE -
NOT FOR DIRECT INFUSION**

Protect from light.

DESCRIPTION

TrophAmine® (10% Amino Acid Injection) is a sterile, nonpyrogenic, hypertonic solution containing crystalline amino acids.

All amino acids designated USP are the “L”-isomer with the exception of Glycine USP, which does not have an isomer.

Each 100 mL contains:

<i>Essential Amino Acids</i>	10%
Isoleucine USP	0.82 g
Leucine USP	1.4 g
Lysine	0.82 g
(added as Lysine Acetate USP)	1.2 g)
Methionine USP	0.34 g
Phenylalanine USP	0.48 g
Threonine USP	0.42 g
Tryptophan USP	0.2 g
Valine USP	0.78 g
Cysteine	<0.016 g
(as Cysteine HCl•H ₂ O USP)	<0.024 g)
Histidine USP*	0.48 g
Tyrosine*	0.24 g
(added as Tyrosine USP	0.044 g
and N-Acetyl-L-Tyrosine	0.24 g)
<i>Nonessential Amino Acids</i>	
Alanine USP	0.54 g
Arginine USP	1.2 g
Proline USP	0.68 g
Serine USP	0.38 g

Glycine USP	0.36 g
L-Aspartic Acid	0.32 g
L-Glutamic Acid	0.5 g
Taurine ^{†‡}	0.025 g
Water for Injection USP	qs
pH adjusted with Glacial Acetic Acid USP	
pH: 5.5 (5.0–6.0)	
Calc. Osmolarity (mOsmol/liter)	866
Total Amino Acids (grams/liter)	100
Total Nitrogen (grams/liter)	15.5
Protein Equivalent (grams/liter)	97

Electrolytes (mEq/liter)

[§] Acetate (CH ₃ COO ⁻)	96.2
Chloride	<3

* Holt LE, Snyderman SE: The amino acid requirements of infants. *JAMA* 1961; 175(2):124-127.

† Rigo J, Senterre J: Is taurine essential for the neonates? *Biol Neonate* 1977; 32:73-76.

‡ Gaul G, Sturman JA, Rähä NCR: Development of mammalian sulfur metabolism: Absence of cystothionase in human fetal tissues. *Pediatr Res* 1972; 6:538-547.

§ Provided as acetic acid and lysine acetate.

CLINICAL PHARMACOLOGY

TrophAmine® provides a mixture of essential and nonessential amino acids as well as taurine and a soluble form of tyrosine, N-Acetyl-L-Tyrosine (NAT). This amino acid composition has been specifically formulated to provide a well tolerated nitrogen source for nutritional support and therapy for infants and young pediatric patients. When administered in conjunction with cysteine hydrochloride, TrophAmine® results in the normalization of the plasma amino acid concentrations to a profile consistent with that of a breast-fed infant.

The rationale for TrophAmine® (Amino Acid Injection) is based on the observation of inadequate levels of essential amino acids in the plasma of infants receiving total parenteral nutrition (TPN) using conventional amino acid solutions. The TrophAmine® formula was developed through the application of specific pharmacokinetic multiple regression analysis relating amino acid intake to the resulting plasma amino acid concentrations.

Clinical studies in infants and young pediatric patients who required TPN therapy showed that infusion of TrophAmine® with a cysteine hydrochloride admixture resulted in a normalization of the plasma amino acid concentrations. In addition, weight gains, nitrogen balance, and serum protein concentrations were consistent with an improving nutritional status.

When infused with hypertonic dextrose as a calorie source, supplemented with cysteine hydrochloride, electrolytes, vitamins, and minerals, TrophAmine® provides total parenteral nutrition in infants and young pediatric patients, with the exception of essential fatty acids.

It is thought that the acetate from lysine acetate and acetic acid, under the conditions of parenteral nutrition, does not impact net acid-base balance when renal and respiratory functions are normal. Clinical evidence seems to support this thinking; however, confirmatory experimental evidence is not available.

The amount of chloride present in TrophAmine® is not of clinical significance.

The addition of cysteine hydrochloride will contribute to the chloride load.

The electrolyte content of any additives that are introduced should be carefully considered and included in total input computations.

INDICATIONS AND USAGE

TrophAmine® is indicated for the nutritional support of infants (including those of low birth weight) and young pediatric patients requiring TPN via either central or peripheral infusion routes. Parenteral nutrition with TrophAmine® is indicated to prevent nitrogen and weight loss or treat negative nitrogen balance in infants and young pediatric patients where (1) the alimentary tract, by the oral, gastrostomy, or jejunostomy route, cannot or should not be used, or adequate protein intake is not feasible by these routes; (2) gastrointestinal absorption of protein is impaired; or (3) protein requirements are substantially increased as with extensive burns. Dosage, route of administration, and concomitant infusion of non-protein calories are dependent on various factors, such as nutritional and metabolic status of the patient, anticipated duration of parenteral nutritional support, and vein tolerance. See **WARNINGS, PRECAUTIONS, Pediatric Use, AND DOSAGE AND ADMINISTRATION.**

Central Venous Nutrition

Central venous infusion should be considered when amino acid solutions are to be admixed with hypertonic dextrose to promote protein synthesis in hypercatabolic or severely depleted infants, or those requiring long-term parenteral nutrition.

Peripheral Parenteral Nutrition

For moderately catabolic or depleted patients in whom the central venous route is not indicated, diluted amino acid solutions mixed with 5-10% dextrose solutions may be infused by peripheral vein, supplemented, if desired, with fat emulsion. In pediatric patients, the final solution should not exceed twice normal serum osmolarity (718 mOsmol/L).

CONTRAINDICATIONS

TrophAmine® is contraindicated in patients with untreated anuria, hepatic coma, inborn errors of amino acid metabolism, including those involving branched chain amino acid metabolism such as maple syrup urine disease and isovaleric acidemia, or hypersensitivity to one or more amino acids present in the solution.

WARNINGS

Safe, effective use of parenteral nutrition requires a knowledge of nutrition as well as

clinical expertise in recognition and treatment of the complications which can occur. **Frequent clinical evaluation and laboratory determinations are necessary for proper monitoring of parenteral nutrition.** Studies should include blood sugar, serum proteins, kidney and liver function tests, electrolytes, hemogram, carbon dioxide content, serum osmolalities, blood cultures, and blood ammonia levels.

WARNING: This product contains aluminum that may be toxic. Aluminum may reach toxic levels with prolonged parenteral administration if kidney function is impaired. Premature neonates are particularly at risk because their kidneys are immature, and they require large amounts of calcium and phosphate solutions, which contain aluminum.

Research indicates that patients with impaired kidney function, including premature neonates, who receive parenteral levels of aluminum at greater than 4 to 5 mcg/kg/day accumulate aluminum at levels associated with central nervous system and bone toxicity. Tissue loading may occur at even lower rates of administration.

Administration of amino acids in the presence of impaired renal function or gastrointestinal bleeding may augment an already elevated blood urea nitrogen. Patients with azotemia from any cause should not be infused with amino acids without regard to total nitrogen intake.

Administration of intravenous solutions can cause fluid and/or solute overload resulting in dilution of serum electrolyte concentrations, overhydration, congested states, or pulmonary edema. The risk of dilutional states is inversely proportional to the electrolyte concentrations of the solutions. The risk of solute overload causing congested states with peripheral and pulmonary edema is directly proportional to the electrolyte concentrations of the solutions.

Administration of amino acid solutions to a patient with hepatic insufficiency may result in plasma amino acid imbalances, hyperammonemia, prerenal azotemia, stupor and coma.

Hyperammonemia is of **special significance in infants** as its occurrence in the syndrome caused by genetic metabolic defects is sometimes associated, although not necessarily in a causal relationship, with mental retardation. This reaction appears to be dose related and is more likely to develop during prolonged therapy. It is essential that blood ammonia be measured frequently in infants. The mechanisms of this reaction are not clearly defined but may involve genetic defects and immature or subclinically impaired liver function.

Conservative doses of amino acids should be given, dictated by the nutritional status of the patient. Should symptoms of hyperammonemia develop, amino acid administration should be discontinued and the patient's clinical status reevaluated.

PRECAUTIONS

General

Clinical evaluation and periodic laboratory determinations are necessary to monitor changes in fluid balance, electrolyte concentrations, and acid-base balance during prolonged parenteral therapy or whenever the condition of the patient warrants such evaluation. Significant deviations from normal concentrations may require the use of

additional electrolyte supplements.

Strongly hypertonic nutrient solutions should be administered via an intravenous catheter placed in a central vein, preferably the superior vena cava.

Care should be taken to avoid circulatory overload, particularly in patients with cardiac insufficiency.

Special care must be taken when giving hypertonic dextrose to a diabetic or pre-diabetic patient. To prevent severe hyperglycemia in such patients, insulin may be required.

Administration of glucose at a rate exceeding the patient's utilization rate may lead to hyperglycemia, coma, and death.

Administration of amino acids without carbohydrates may result in the accumulation of ketone bodies in the blood. Correction of this ketonemia may be achieved by the administration of carbohydrate.

Peripheral administration of TrophAmine® (Amino Acid Injection) requires appropriate dilution and provision of adequate calories. Care should be taken to assure proper placement of the needle within the lumen of the vein. The venipuncture site should be inspected frequently for signs of infiltration. If venous thrombosis or phlebitis occurs, discontinue infusions or change infusion site and initiate appropriate treatment. In pediatric patients, the final solution should not exceed twice normal serum osmolarity (718 mOsm/L).

Extraordinary electrolyte losses such as may occur during protracted nasogastric suction, vomiting, diarrhea, or gastrointestinal fistula drainage may necessitate additional electrolyte supplementation.

Metabolic acidosis can be prevented or readily controlled by adding a portion of the cations in the electrolyte mixture as acetate salts and in the case of hyperchloremic acidosis, by keeping the total chloride content of the infusate to a minimum.

TrophAmine® (Amino Acid Injection) contains less than 3 mEq chloride per liter.

TrophAmine® contains no added phosphorus. Patients, especially those with hypophosphatemia, may require the addition of phosphate. To prevent hypocalcemia, calcium supplementation should always accompany phosphate administration. To assure adequate intake, serum levels should be monitored frequently.

To minimize the risk of possible incompatibilities arising from mixing this solution with other additives that may be prescribed, the final infusate should be inspected for cloudiness or precipitation immediately after mixing, prior to administration, and periodically during administration.

Use only if solution is clear, the seal unbroken, and vacuum is present.

Drug product contains no more than 25 mcg/L of aluminum.

Laboratory Tests

Frequent clinical evaluation and laboratory determinations are necessary for proper monitoring during administration.

Laboratory tests should include measurement of blood sugar, electrolyte, and serum protein concentrations; kidney and liver function tests; and evaluation of acid-base balance and fluid balance. Other laboratory tests may be suggested by the patient's

condition.

Drug Interactions

Some additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic techniques. Mix thoroughly. Do not store.

Carcinogenesis, Mutagenesis, Impairment of Fertility

No *in vitro* or *in vivo* carcinogenesis, mutagenesis, or fertility studies have been conducted with TrophAmine®.

Pregnancy

Teratogenic Effects

Animal reproduction studies have not been conducted with TrophAmine® (Amino Acid Injection). It is also not known whether TrophAmine® can cause fetal harm when administered to a pregnant woman or can affect reproduction capacity. TrophAmine® should be given to a pregnant woman only if clearly needed.

Labor and Delivery

Information is unknown.

Nursing Mothers

It is not known whether this drug is excreted in human milk. Because many drugs are excreted in human milk, caution should be exercised with TrophAmine® if administered to a nursing woman.

Pediatric Use

As in all cases of fluid and electrolyte replacement and parenteral nutrition, careful monitoring and special caution is required in pediatric use, especially in pediatric patients with renal failure, acute sepsis, or low birth weight.

The total volume of nutritional fluid and the rate of administration in each patient will be based on individually calculated maintenance and/or replacement fluid requirements, and nutritional needs, and will vary with the child's age, body weight and renal function.

In neonates and very small infants, particularly careful monitoring will be required to maintain fluid and electrolyte balance, including monitoring of blood glucose.

See **INDICATIONS AND USAGE, WARNINGS, and DOSAGE AND ADMINISTRATION.**

Geriatric Use

TrophAmine® has not been studied in geriatric patients. Elderly patients are known to be more prone to fluid overload and electrolyte imbalance than younger patients. This may be related to impairment of renal function which is more frequent in an elderly population. As a result the need for careful monitoring of fluid and electrolyte therapy is greater in the elderly.

All patients, including the elderly, require an individual dose of all parenteral nutrition products to be determined by the physician on an individual case-by-case basis, which will be based on body weight, clinical condition and the results of laboratory monitoring tests. There is no specific geriatric dose.

See **WARNINGS**.

Special Precautions for Central Venous Nutrition

Administration by central venous catheter should be used only by those familiar with this technique and its complications.

Central venous nutrition may be associated with complications which can be prevented or minimized by careful attention to all aspects of the procedure, including solution preparation, administration, and patient monitoring. **It is essential that a carefully prepared protocol, based on current medical practices, be followed, preferably by an experienced team.**

Although a detailed discussion of the complications is beyond the scope of this insert, the following summary lists those based on current literature:

Technical. The placement of a central venous catheter should be regarded as a surgical procedure. One should be fully acquainted with various techniques of catheter insertion as well as recognition and treatment of complications. For details of techniques and placement sites, consult the medical literature. X-ray is the best means of verifying catheter placement. Complications known to occur from the placement of central venous catheters are pneumothorax, hemothorax, hydrothorax, artery puncture and transection, injury to the brachial plexus, malposition of the catheter, formation of arteriovenous fistula, phlebitis, thrombosis, and air and catheter embolus.

Septic. The constant risk of sepsis is present during central venous nutrition. Since contaminated solutions and infusion catheters are potential sources of infection, it is imperative that the preparation of parenteral nutrition solutions and the placement and care of catheters be accomplished under controlled aseptic conditions.

Solutions should ideally be prepared in the hospital pharmacy in a laminar flow hood. The key factor in their preparation is careful aseptic technique to avoid inadvertent touch contamination during mixing of solutions and subsequent admixtures.

Parenteral nutrition solutions should be used promptly after mixing. Any storage should be under refrigeration for as brief a time as possible. Administration time for a TPN solution should not exceed 24 hours.

Consult the medical literature for a discussion of the management of sepsis during central venous nutrition. In brief, typical management includes replacing the solution being administered with a fresh container and set, and the remaining contents are cultured for bacterial or fungal contamination. If sepsis persists and another source of infection is not identified, the catheter is removed, the proximal tip cultured, and a new catheter reinserted when the fever has subsided. Non-specific, prophylactic antibiotic treatment is not recommended. Clinical experience indicates that the catheter is likely to be the prime source of infection as opposed to aseptically prepared and properly stored solutions.

Metabolic. The following metabolic complications have been reported: metabolic acidosis, hypophosphatemia, alkalosis, hyperglycemia and glycosuria, osmotic diuresis

and dehydration, rebound hypoglycemia, elevated liver enzymes, hypo- and hypervitaminosis, electrolyte imbalances, and hyperammonemia in pediatric patients. Frequent clinical evaluation and laboratory determinations are necessary, especially during the first few days of venous nutrition, to prevent or minimize these complications.

ADVERSE REACTIONS

See “**WARNINGS**” and “***Special Precautions for Central Venous Nutrition.***”

Reactions reported in clinical studies as a result of infusion of the parenteral fluid were water weight gain, edema, increase in BUN, and mild acidosis.

Reactions which may occur because of the solution or the technique of administration include febrile response, infection at the site of injection, venous thrombosis or phlebitis extending from the site of injection, extravasation and hypervolemia.

Local reaction at the infusion site, consisting of a warm sensation, erythema, phlebitis and thrombosis, have been reported with peripheral amino acid infusions, especially if other substances are also administered through the same site.

If electrolyte supplementation is required during peripheral infusion, it is recommended that additives be administered throughout the day in order to avoid possible venous irritation. Irritating additive medications may require injection at another site and should not be added directly to the amino acid infusate.

Symptoms may result from an excess or deficit of one or more of the ions present in the solution; therefore, frequent monitoring of electrolyte levels is essential.

Phosphorus deficiency may lead to impaired tissue oxygenation and acute hemolytic anemia. Relative to calcium, excessive phosphorus intake can precipitate hypocalcemia with cramps, tetany and muscular hyperexcitability.

If an adverse reaction does occur, discontinue the infusion, evaluate the patient, institute appropriate therapeutic countermeasures and save the remainder of the fluid for examination if deemed necessary.

OVERDOSAGE

In the event of a fluid or solute overload during parenteral therapy, reevaluate the patient's condition, and institute appropriate corrective treatment.

DOSAGE AND ADMINISTRATION

The objective of nutritional management of infants and young pediatric patients is the provision of sufficient amino acid and caloric support for protein synthesis and growth.

The total daily dose of TrophAmine® (Amino Acid Injection) depends on daily protein requirements and on the patient's metabolic and clinical response. The determination of nitrogen balance and accurate daily body weights, corrected for fluid balance, are probably the best means of assessing individual protein requirements. Dosage should also be guided by the patient's fluid intake limits and glucose and nitrogen tolerances, as well as by metabolic and clinical response.

Recommendations for allowances of protein in infant nutrition have ranged from 2 to 4 grams of protein per kilogram of body weight per day (2 to 4 g/kg/day).¹The recommended dosage of TrophAmine® is 2 to 2.5 grams of amino acids per kilogram of body weight per day (2 to 2.5 g/kg/day) for infants up to 10 kilograms. For infants and young pediatric patients larger than 10 kilograms, the total dosage of amino acids should include the 20 to 25 grams/day for the first 10 kg of body weight plus 1 to 1.25 g/day for each kg of body weight over 10 kilograms.

Typically, TrophAmine® is admixed with 50% or 70% Dextrose Injection USP supplemented with electrolytes and vitamins and administered over a period of time not to exceed 24 hours.

Total daily fluid intake should be appropriate for the patient's age and size. A fluid dose of 125 mL per kilogram body weight per day is appropriate for most infants on TPN. Although nitrogen requirements may be higher in severely hypercatabolic or depleted patients, provision of additional nitrogen may not be possible due to fluid intake limits, nitrogen, or glucose intolerance.

Cysteine is considered to be an essential amino acid in infants and young pediatric patients. An admixture of cysteine hydrochloride to the TPN solution is therefore recommended. Based on clinical studies, the recommended dosage is 1 mmole of L-cysteine hydrochloride monohydrate per kilogram of body weight per day.

In many patients, provision of adequate calories in the form of hypertonic dextrose may require the administration of exogenous insulin to prevent hyperglycemia and glycosuria. To prevent rebound hypoglycemia, a solution containing 5% dextrose should be administered when hypertonic dextrose solutions are abruptly discontinued.

Fat emulsion coadministration should be considered when prolonged (more than 5 days) parenteral nutrition is required in order to prevent essential fatty acid deficiency (E.F.A.D.). Serum lipids should be monitored for evidence of E.F.A.D. in patients maintained on fat free TPN.

The provision of sufficient intracellular electrolytes, principally potassium, magnesium, and phosphate, is required for optimum utilization of amino acids. In addition, sufficient quantities of the major extracellular electrolytes sodium, calcium, and chloride, must be given. In patients with hyperchloremic or other metabolic acidoses, sodium and potassium may be added as the acetate salts to provide bicarbonate precursor. The electrolyte content of TrophAmine® must be considered when calculating daily electrolyte intake. Serum electrolytes, including magnesium and phosphorus, should be monitored frequently.

Appropriate vitamins, minerals and trace elements should also be provided.

Central Venous Nutrition. Hypertonic mixtures of amino acids and dextrose may be safely administered by continuous infusion through a central venous catheter with the tip located in the superior vena cava. Initial infusion rates should be slow, and gradually increased to the recommended 60-125 mL per kilogram body weight per day. If administration rate should fall behind schedule, no attempt to "catch up" to planned intake should be made. In addition to meeting protein needs, the rate of administration, particularly during the first few days of therapy, is governed by the patient's glucose tolerance. Daily intake of amino acids and dextrose should be increased gradually to the maximum required dose as indicated by frequent determinations of glucose levels in

blood and urine.

Use of an infusion pump is advisable to maintain a steady infusion rate during central venous infusion.

Peripheral Parenteral Nutrition. For patients in whom the central venous route is not indicated and who can consume adequate calories enterally, TrophAmine® (Amino Acid Injection) may be administered by peripheral vein with or without parenteral carbohydrate calories. Such infusates can be prepared by dilution with B. Braun's Sterile Water for Injection or 5%-10% Dextrose Injection to prepare isotonic or slightly hypertonic solutions for peripheral infusion. It is essential that peripheral infusion be accompanied by adequate caloric intake. In pediatric patients, the final solution should not exceed twice normal serum osmolarity (718 mOsmol/L).

Parenteral drug products should be inspected visually for particulate matter and discoloration prior to administration, whenever solution and container permit.

TrophAmine® may be admixed with solutions which contain phosphate or which have been supplemented with phosphate. The presence of calcium and magnesium ions in an additive solution should be considered when phosphate is also present, in order to avoid precipitation.

Care must be taken to avoid incompatible admixtures. Consult with pharmacist.

1 Suskind RM: Textbook of Pediatric Nutrition, Raven Press, New York, 1981.

HOW SUPPLIED

TrophAmine® pharmacy bulk package is supplied sterile and nonpyrogenic in 500 mL flexible plastic bags, packaged 12 per case.

NDC	REF	Size
0264-1933-10	S9333	500 mL

Not made with natural rubber latex, DEHP or PVC.

Storage

Store in the closed corrugated case; do not expose TrophAmine® and/or TrophAmine®-containing admixed parenteral nutrition (PN) solution to light. Do not remove container from overwrap until ready to use. Do not use if overwrap has been opened or damaged. Exposure of pharmaceutical products to heat should be minimized. Avoid excessive heat. It is recommended that the product be stored at 20°C to 25°C (68°F to 77°F); excursions permitted between 15°C to 30°C (59°F to 86°F). [See USP Controlled Room Temperature]. Solution that has been frozen must not be used.

Rx only

TrophAmine is a registered trademark of B. Braun Medical Inc.

Preparation Instructions for TrophAmine®


Caution: TrophAmine® Pharmacy Bulk Package is intended for intravenous infusion after further dilution only.

Do not expose TrophAmine® and/or TrophAmine®-containing admixed parenteral nutrition (PN) solution to light.


The pharmacy bulk package is for use in a Pharmacy Admixture Service only. Use of this product is restricted to a suitable work area, such as a laminar flow hood (or an equivalent clean air compounding area).

Designed for use with an unvented sterile dispensing set.


Step 1: Inspect Pharmacy Bulk Package overwrap and primary bag and do not use if damaged. Inspect oxygen indicator and do not use if oxygen indicator is pink or dark pink. Use only if container and seals are intact.




Step 2: To open, tear overwrap starting from the tear notches (Figure 1). Remove TrophAmine® Pharmacy Bulk Package bag from overwrap and discard oxygen indicator, oxygen absorber and overwrap.




Step 3: Inspect TrophAmine® Pharmacy Bulk Package bag visually (Figure 2). Parenteral drug products should be inspected visually for particulate matter and discoloration prior to administration, whenever solution and container permit. Inspect TrophAmine® for cloudiness, haze or particulate matter in good lighting. Discard the bag if any particulates or discoloration are observed.





Step 4: Remove aluminum foil of set port at the bottom of the Pharmacy Bulk Package bag (Figure 3a) and attach an unvented sterile dispensing set (Figure 3b).



The set port may be penetrated only one time, utilizing a suitable sterile dispensing set which allows measured dispensing of the contents.

The withdrawal of container contents should be accomplished without delay using aseptic technique. Discard container no later than 4 hours after initial closure puncture.

Use of a syringe with needle is not recommended. Multiple entries increase the potential of the microbial and particulate contamination.

B. Braun Medical Inc.

Bethlehem, PA 18018-3524 USA

1-800-227-2862

Made in Germany

Revised: June 2022

LD-625-2

PRINCIPAL DISPLAY PANEL - 500 mL Container Label

**TrophAmine®
(10% Amino Acid Injection)**

NDC 0264-1933-10

REF S9333

500 mL

Protect from light.

For Intravenous Infusion After Dilution

852/12627204/0522

LOT

EXP:

SN:

Recommended Storage: Store at 20°C to 25°C (68°F to 77°F). [See USP Controlled Room Temperature]. Solution that has been frozen must not be used. Do not expose to light. See package insert for Dosage and Administration.

PHARMACY BULK PACKAGE - NOT FOR DIRECT INFUSION

Each 100 mL contains:

Essential Amino Acids	
Isoleucine USP	0.82 g
Leucine USP	1.4 g
Lysine	0.82 g
(added as Lysine Acetate USP	1.2 g)
Methionine USP	0.34 g
Phenylalanine USP	0.48 g
Threonine USP	0.42 g
Tryptophan USP	0.2 g
Valine USP	0.78 g
Cysteine	<0.016 g
(as Cysteine HCl•H ₂ O USP	<0.024 g)
Histidine USP	0.48 g
Tyrosine	0.24 g
(added as Tyrosine USP	0.044 g
and N-Acetyl-L-Tyrosine	0.24 g)

Affix accompanying label for date and time of entry

Nonessential Amino Acids

Alanine USP	0.54 g
Arginine USP	1.2 g
Proline USP	0.68 g
Serine USP	0.38 g
Glycine USP	0.36 g
L-Aspartic Acid	0.32 g
L-Glutamic Acid	0.5 g
Taurine	0.025 g

Water for Injection USP qs
pH adjusted with Glacial Acetic Acid USP
pH: 5.5 (5.0-6.0)
Calc. Osmolarity: 866 mOsmol/liter

Electrolytes (mEq/liter):

Chloride <3

Acetate 96.2 (see Package Insert)

Sterile, nonpyrogenic.

Single dose container.

Use only if bag and seal are undamaged and solution is clear.

Do not remove container from overwrap until ready to use.



Do not use if overwrap has been opened or damaged.

For intravenous infusion after dilution.

Not made with natural rubber latex, DEHP or PVC.

Rx only

TrophAmine is a registered trademark of B. Braun Medical Inc.

B. Braun Medical Inc.

Bethlehem, PA 18018-3524 USA

1-800-227-2862

Made in Germany

LD-626-2

TrophAmine[®]

(10% Amino Acid Injection)

500 mL

Protect from light.

For Intravenous Infusion After Dilution
852/12627204/0522

NDC 0264-1933-10

REF S9333

LOT

EXP:

SN:

Recommended Storage:

Store at 20°C to 25°C (68°F to 77°F).

[See USP Controlled Room Temperature].

Solution that has been frozen must not be used. Do not expose to light. See package insert for Dosage and Administration.

400 mL

**PHARMACY BULK PACKAGE -
NOT FOR DIRECT INFUSION**

Affix accompanying label
for date and time of entry

300 mL

Each 100 mL contains:

Essential Amino Acids

Isoleucine USP	0.82 g
Leucine USP	1.4 g
Lysine	0.82 g
(added as Lysine Acetate USP	1.2 g)
Methionine USP	0.34 g
Phenylalanine USP	0.48 g
Threonine USP	0.42 g
Tryptophan USP	0.2 g
Valine USP	0.78 g

200 mL

Cysteine	<0.016 g
(as Cysteine HCl•H ₂ O USP	<0.024 g)
Histidine USP	0.48 g
Tyrosine	0.24 g
(added as Tyrosine USP	0.044 g
and N-Acetyl-L-Tyrosine	0.24 g)

Nonessential Amino Acids

Alanine USP	0.54 g	Glycine USP	0.36 g
Arginine USP	1.2 g	L-Aspartic Acid	0.32 g
Proline USP	0.68 g	L-Glutamic Acid	0.5 g
Serine USP	0.38 g	Taurine	0.025 g

Water for Injection USP qs

pH adjusted with Glacial Acetic Acid USP

pH: 5.5 (5.0-6.0)

Calc. Osmolarity: 866 mOsmol/liter

Electrolytes (mEq/liter):

Chloride <3

Acetate 96.2 (see Package Insert)

Sterile, nonpyrogenic.

Single dose container.

Use only if bag and seal are undamaged and solution is clear.

Do not remove container from overwrap until ready to use.



B. Braun Medical Inc.
Bethlehem, PA 18018-3524 USA
1-800-227-2862
Made in Germany

LD-626-2

Do not use if overwrap has been opened or damaged.

For intravenous infusion after dilution.

Not made with natural rubber latex, DEHP or PVC. Rx only

TrophAmine is a registered trademark of B. Braun Medical Inc.

PRINCIPAL DISPLAY PANEL - Accompanying 500 mL Bag Label

Date of Entry: __/__/____

Time of Entry: _____ AM/PM

Discard _____ hours after initial entry

LD-628-2

852/12627207/0522

Date of Entry: __/__/____

Time of Entry: _____ AM/PM

Discard _____ hours after initial entry

LD-628-2
852/12627207/0522

TROPHAMINE

isoleucine, leucine, lysine acetate, methionine, phenylalanine, threonine, tryptophan, valine, cysteine hydrochloride, histidine, tyrosine, n-acetyl-tyrosine, alanine, arginine, proline, serine, glycine, aspartic acid, glutamic acid, and taurine solution

Product Information

Product Type	HUMAN PRESCRIPTION DRUG	Item Code (Source)	NDC:0264-1933
Route of Administration	INTRAVENOUS		

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
ISOLEUCINE (UNII: 04Y7590D77) (ISOLEUCINE - UNII:04Y7590D77)	ISOLEUCINE	0.82 g in 100 mL
LEUCINE (UNII: GMW67QNF9C) (LEUCINE - UNII:GMW67QNF9C)	LEUCINE	1.4 g in 100 mL
LYSINE ACETATE (UNII: TTL6G7LIWZ) (LYSINE - UNII:K3Z4F929H6)	LYSINE	1.2 g in 100 mL
METHIONINE (UNII: AE28F7PNPL) (METHIONINE - UNII:AE28F7PNPL)	METHIONINE	0.34 g in 100 mL

PHENYLALANINE (UNII: 47E5O17Y3R) (PHENYLALANINE - UNII:47E5O17Y3R)	PHENYLALANINE	0.48 g in 100 mL
THREONINE (UNII: 2ZD004190S) (THREONINE - UNII:2ZD004190S)	THREONINE	0.42 g in 100 mL
TRYPTOPHAN (UNII: 8DUH1N11BX) (TRYPTOPHAN - UNII:8DUH1N11BX)	TRYPTOPHAN	0.2 g in 100 mL
VALINE (UNII: HG18B9YRS7) (VALINE - UNII:HG18B9YRS7)	VALINE	0.78 g in 100 mL
CYSTEINE HYDROCHLORIDE (UNII: ZT934N0X4W) (CYSTEINE - UNII:K848JZ4886)	CYSTEINE	0.024 g in 100 mL
HISTIDINE (UNII: 4QD397987E) (HISTIDINE - UNII:4QD397987E)	HISTIDINE	0.48 g in 100 mL
TYROSINE (UNII: 42HK56048U) (TYROSINE - UNII:42HK56048U)	TYROSINE	0.044 g in 100 mL
ACETYL L-TYROSINE (UNII: DA8G610ZO5) (ACETYL L-TYROSINE - UNII:DA8G610ZO5)	ACETYL L-TYROSINE	0.24 g in 100 mL
ALANINE (UNII: OF5P57N2ZX) (ALANINE - UNII:OF5P57N2ZX)	ALANINE	0.54 g in 100 mL
ARGININE (UNII: 94ZLA3W45F) (ARGININE - UNII:94ZLA3W45F)	ARGININE	1.2 g in 100 mL
PROLINE (UNII: 9DLQ4CIU6V) (PROLINE - UNII:9DLQ4CIU6V)	PROLINE	0.68 g in 100 mL
SERINE (UNII: 452VLY9402) (SERINE - UNII:452VLY9402)	SERINE	0.38 g in 100 mL
GLYCINE (UNII: TE7660XO1C) (GLYCINE - UNII:TE7660XO1C)	GLYCINE	0.36 g in 100 mL
ASPARTIC ACID (UNII: 30KYC7MIAI) (ASPARTIC ACID - UNII:30KYC7MIAI)	ASPARTIC ACID	0.32 g in 100 mL
GLUTAMIC ACID (UNII: 3KX376GY7L) (GLUTAMIC ACID - UNII:3KX376GY7L)	GLUTAMIC ACID	0.5 g in 100 mL
TAURINE (UNII: 1EQV5MLY3D) (TAURINE - UNII:1EQV5MLY3D)	TAURINE	0.025 g in 100 mL

Inactive Ingredients

Ingredient Name	Strength
WATER (UNII: 059QF0K00R)	
ACETIC ACID (UNII: Q40Q9N063P)	

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:0264-1933-10	12 in 1 CASE	10/27/2020	
1		500 mL in 1 CONTAINER; Type 0: Not a Combination Product		

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
NDA	NDA019018	10/27/2020	

Labeler - B. Braun Medical Inc. (002397347)

